

Medcare Lift Slings



Care Sling

This is our most popular sling style and is designed to fit the needs of most patients. It is available in ripstop polyester for regular use, and mesh for bathing.

Polyester:

- 400018: Small Sling
- 400020: Medium Sling
- 400022: Large Sling
- 400017: X-Large Sling
- 400015: Bariatric Sling

Mesh:

- 400008: Small Sling
- 400010: Medium Sling
- 400012: Large Sling
- 400011: X-Large Sling
- 400014: Bariatric Sling



Our wide range of sling styles meets many multicare needs. Medcare slings slightly stretch to provide a contoured fit to most body shapes, while the 4-point hook-up system ensures patient comfort.

Our polyester slings are made from our proven ripstop polyester material, which is rated and tested to 1,000 lbs. and if punctured or torn, it will not rip!

Care Sling with Head Support

This sling is ideal for patients with weak or poor head control. Head support slings are available in polyester for regular use, or mesh for bathing.

Polyester:


- 400028: Small Sling w/ Head Support
- 400030: Medium Sling w/ Head Support
- 400032: Large Sling w/ Head Support
- 400034: X-Large Sling w/ Head Support
- 400036: Bariatric Sling w/ Head Support

Mesh:


- 400038: Small Sling w/ Head Support
- 400040: Medium Sling w/ Head Support
- 400042: Large Sling w/ Head Support
- 400043: X-Large Sling w/ Head Support
- 400041: Bariatric Sling w/ Head Support



Sling Fabric Swatches



Ripstop Polyester



Mesh

Medcare Slings Contain Antimicrobial Protection:

Antimicrobial protection works 24/7 for an added level of cleanliness protection and inhibits growth of microbes such as bacteria, mold, and mildew that can cause stains, odors, and product deterioration on slings.

Repositioning Sling

Designed to be left on the bed under the sheet, this breathable sling, along with a Medcare Ceiling Lift, makes lateral transfer and turning or boosting a patient very easy.

Polyester:

- 402503: Bariatric Sling
Dimensions: 70" W x 88" L

Mesh:

- 402500: Repositioning Sling
Dimensions: 55.5" W x 80.5" L



Commode Sling

Designed to provide easy transfers to and from commodes, toilets, or for hygienic purposes. An enlarged opening makes utilizing the commode quicker and easier than standard slings with minimal interference of clothing or sling material. An adjustable belt attaches around the waist for added security and comfort.

Mesh:

- 400071: Small Commode Sling
- 400072: Medium Commode Sling
- 400073: Large Commode Sling
- 400074: X-Large Commode Sling



Stretcher Sling

For patients that must be kept in a supine position, this sling allows easy transfer of a patient from bed to gurney or operating table. It is breathable and may be left under patient.

Polyester:

- 402450: Stretcher Sling
Dimensions: 47.5" W x 80.5" L



Lift Slings (cont.)

Multi-Purpose Sling

Designed for a variety of uses, Multi-Purpose Slings have wider straps that are more comfortable between the legs. These straps can also be configured so straps do not go between legs.

Polyester:

- 400024: Small Multi-Purpose Sling
- 400025: Medium Multi-Purpose Sling
- 400026: Large Multi-Purpose Sling
- 400027: X-Large Multi-Purpose Sling

Mesh:

- 400044: Small Multi-Purpose Sling
- 400046: Medium Multi-Purpose Sling
- 400048: Large Multi-Purpose Sling
- 400047: X-Large Multi-Purpose Sling



Sling Sizing:


Use the sizing chart below as a general guide. Keep in mind that two patients that are the same weight may have a completely different body shape and size and may require different size slings.

Size/Strap Color	Weight
Small (Black) ●	0–100 lbs.
Medium (Blue) ●	100–210 lbs.
Large (Yellow) ●	210–310 lbs.
X-Large (White) ○	310–450 lbs.
Bariatric (Purple) ●	450–1000 lbs.



Washing Instructions:

- Machine wash in soap solution
- Max washing temperature: 165° F
- Rinse thoroughly
- Do not bleach
- Do not dry clean



Drying Instructions:

- Hang or tumble dry
- Max drying temperature: 175° F
- Do not overload dryer
- Do not place with direct heat source in dryer

Limb Sling

Designed to assist in lifting a leg or arm for cleaning, changing wound dressings, catheter insertion, etc. Available in two widths: 6" and 12".

Polyester:

- 400250: Limb Sling 6"
- 400275: Limb Sling 12"



Long Seated Sling

Designed to lift a patient much like the Care Sling, but without the need to cross the legs. This works great for people with spinal chord injuries, in which a patient needs to be transported in a variety of ways like supine or sitting, etc.

Mesh:

- 400122: Large Long Seated Sling
- 400123: X-Large Long Seated Sling



Medcare Products

151 East Cliff Road
Burnsville, MN 55337

Toll-free: 800 695 4479
Fax: 952 300 2100

www.MedcareLifts.com

© 2014 Medcare Products Inc.

Rev 031814

IMPORTANT: Always assess the patient to ensure the right fit and comfort level. Sling type used for different lifting situations is always at the discretion of the facility. For tips and training on using Medcare's various slings, contact your local Safe Patient Handling Consultant or Medcare's technical support at 1-800-695-4479.