

All slings are rated to 1,000 lbs. capacity. Sizing recommendations are provided for patient comfort and fit.

<u>Patient Weight</u>	<u>Size/Strap Color</u>
0 - 100 lbs. 0 - 45 kg	Small Black
100 - 210 lbs. 45 - 95 kg	Medium Blue
210 - 310 lbs. 95 - 141 kg	Large Yellow
310 - 450 lbs. 141 - 204 kg	X-Large White
450 - 1000 lbs. 295 - 454 kg	Bariatric Purple



Photo of Medium Care Sling without Head Support (400020)

Disposable Care Sling

All slings are rated to 1,000 lbs. capacity. Sizing recommendations are provided for patient comfort and fit.

Patient Weight Size/Strap Color

0 - 100 lbs. 0 - 45 kg	Small	Black
100 - 210 lbs. 45 - 95 kg	Medium	Blue
210 - 310 lbs. 95 - 141 kg	Large	Yellow
310 - 450 lbs. 141 - 204 kg	X-Large	White
450 - 1000 lbs. 295 - 454 kg	Bariatric	Purple



Photo of Medium Care Sling without Head Support (600020)