

# WE'RE THE COMPANY THAT BRINGS YOU The Independent Lifter

TRANSFER AID



The Independent Lifter was designed to provide independence and freedom to various transfer situations. In combination with a Handicare ceiling lift, the Independent Lifter removes the need for a sling, and gives simplicity and accessibility for convenient application for situations including swimming pools, toilets, showers, and others. The Independent Lifter is perfect for homecare where simplicity, independence, and safety are needed.

## SPECIFICATIONS

- Width: 20" (Adjustable)
- Height: 30"
- Length: 23"
- Weight limit: 440 lbs

## PRODUCT HIGHLIGHTS

- Easily attaches to Handicare ceiling lifts with a QRS Hook
- Padding to provide support and comfort for the user
- Stainless steel design for longevity and reliability
- Adjustable leg straps for maximum comfort and preferred positioning for patient
- Three position spring mechanism to adjust gripping strength of lifter for different body types
- Leg Supports come in three different sizes: Small, Medium, Large

