

# Care Stand Belts & Accessories

Medcare's Care Stand Belts are for weight-bearing patients/residents who have more consistent and reliable motor functions. Care Stand Belts are the first and only sit-to-stand belts available with Microban antimicrobial product protection. Designed for use with Care Stand, Care Bar Stand, Care Stand-n-Weigh and Care Bar Stand-n-Weigh.

## Medcare Disposable Slings & Belts Contain Antimicrobial Product Protection

Microban protection provides continuous antimicrobial cleaning that inhibits the growth of microbes such as bacteria, mold and mildew that can cause stains, odors and product deterioration, keeping the slings cleaner between cleanings.





## Care Belt

Care belts are designed to go behind the resident's back and buckle around the waist. Reinforced, color-coded straps allow for easy adjustment and safety. These belts are for weight-bearing residents who have more consistent and reliable motor functions.

### Item Numbers:

- 400078 – Small Care Belt: Fits waists 24" to 48"
- 400080 – Medium Care Belt: Fits waists 30" to 54"
- 400082 – Large Care Belt: Fits waists 36" to 60"
- 400084 – X-Large Care Belt: Fits waists 42" to 66"



## Sheepskin Cover for Care Stand Belt

Used in conjunction with the Stand Belt, the sheepskin cover provides an added level of comfort and reduces slippage, keeping the belt in place.

### Item Numbers:

- 400106 – Small Sheepskin Cover
- 400100 – Medium Sheepskin Cover
- 400102 – Large Sheepskin Cover
- 400104 – X-Large Sheepskin Cover
- 400085 – Belt & Cover Combo
- 400086 – Large Stand Belt and Sheepskin Cover Combination



## Mesh Sling Bag

Attaches to your Medcare Lift or Medcare Stand for convenient storage of your slings and belts.

### Item Number:

- 400094 – Mesh Sling Bag



## Medcare Transport Attachment

The Transport Attachment is designed for use in conjunction with the Stand Belt. It provides additional seat and back support when needed. The Transport Attachment can also be used for added support when using the stand to ambulate.



### Item Numbers:

- 400077 – Small Transport Attachment
- 400081 – Medium Transport Attachment
- 400083 – Large Transport Attachment
- 400079 – X-Large Transport Attachment
- 400090 – Medium Belt & Attachment Combo
- 400091 – Large Belt & Attachment Combo



## Belt Sizing Chart

Use this sizing chart as a general guide. Keep in mind that patents/residents who are the same weight may have a completely different body shape and size and may require different size belts.

Belt Size	Belt Waist Size	Strap Color
Small	24" to 48"	Black
Medium	30" to 54"	Blue
Large	36" to 60"	Yellow
X-Large	42" to 66"	White