

Medcare Patient Specific Disposable Slings & Belts



Medcare's Patient Specific Disposable Slings and Belts are designed for use by a single patient and help reduce the risk of cross-contamination. Although manufactured to be disposable, they can be used multiple times on the same patient, and all sizes are rated to hold over 1,000 lbs. Once washed, patient specific slings are no longer safe to use with patients.



“ Medcare's patient specific slings can be used with most mobile ceiling lift styles and fit most 2-, 4-, and 6-point spreader bars. ”




Features:

- Slings hold over 1,000 lbs.
- 100% breathable and may be left under patient
- May be used for bathing, showering, or in the pool as slings will not rip or fall apart when wet
- Can be spot-cleaned or spray disinfected





Sling and Belt Sizing:

Use the sizing chart below as a general guide. Keep in mind that two patients that are the same weight may have a completely different body shape and size and may require different size slings or belts.

Sling Sizes:

Size/Strap Color	Weight
Small (Black) 	0-100 lbs.
Medium (Blue) 	100-210 lbs.
Large (Yellow) 	210-300 lbs.
X-Large (White) 	310-450 lbs.
Bariatric (Purple) 	450-1000 lbs.

Belt Sizes:

Size/Strap Color	Waist Size
Small (Black) 	24 - 48 in.
Medium (Blue) 	30 - 54 in.
Large (Yellow) 	36 - 60 in.
X-Large (White) 	42 - 66 in.



Patient Specific Disposable Slings & Belts (cont.)

Care Sling

This is our most popular sling style and is designed to fit the needs of most patients.

Item Number:

- 600018: Small Sling
- 600020: Medium Sling
- 600022: Large Sling
- 600017: X-Large Sling
- 600015: Bariatric Sling



Disposable Care Sling with Head Support

This sling is ideal for patients with weak or poor head control.

Item Number:

- 600028: Small Sling w/ Head Support
- 600030: Medium Sling w/ Head Support
- 600032: Large Sling w/ Head Support
- 600034: X-Large Sling w/ Head Support
- 600036: Bariatric Sling w/ Head Support



Disposable Limb Sling

Designed to assist in lifting a leg or arm for cleaning, changing wound dressings, catheter insertion, etc. Available in two widths: 6" and 12".

Item Number:

- 600250: Limb Sling 6"
- 600275: Limb Sling 12"



Disposable Commode Sling

Designed to provide easy transfers to and from commodes, toilets, or for hygienic purposes. An enlarged opening makes utilizing the commode quicker and easier than standard slings with minimal interference of clothing or sling material. An adjustable belt attaches around the waist for added security and comfort.

Item Number:

- 600071: Small Commode Sling
- 600072: Medium Commode Sling
- 600073: Large Commode Sling
- 600074: X-Large Commode Sling



Disposable Stand Belt

The Medicare Stand Belt is designed to go behind a patient's back and buckle around the waist. Reinforced, color-coded straps allow for easy adjustment and safety. The belt is for weight-bearing patients who have more consistent and reliable motor functions.

Item Number:

- 600078: Small Stand Belt
- 600080: Medium Stand Belt
- 600082: Large Stand Belt
- 600084: X-Large Stand Belt



Disposable Stretcher Sling

For patients that must be kept in a supine position, this sling allows easy transfer of a patient from bed to gurney or operating table. It is breathable and may be left under patient.

Item Number:

- 600200: Stretcher Sling
- Dimensions: 47.5" W x 80.5" L



Disposable Repositioning Sling

Designed to be left on the bed under the sheet, this breathable sling, along with Medicare Ceiling Lift, makes lateral transfer and turning or boosting a patient very easy.

Item Number:

- 600210: Repositioning Sling
- Dimensions: 55.5" W x 80.5" L



Medcare Products

151 East Cliff Road
Burnsville, MN 55337

Toll-free: 800 695 4479
Fax: 952 300 2100

www.MedcareLifts.com

© 2014 Medcare Products Inc.

Rev 031914

IMPORTANT: Always assess the patient to ensure the right fit and comfort level. Sling type used for different transfers is always at the discretion of the facility. For tips and training on using Medicare's slings and belts, contact your local Safe Patient Handling Consultant or Medicare's technical support at 1-800-695-4479.