



# CARLO Comfort Floor Lift

Safe and Effortless Patient Transfers



## Providing a High Level of **Safety and Comfort**

Lightweight and easy to maneuver the CARLO Comfort is suitable for all areas of care, in addition to lifting and transporting the client the lift can also be used for bathing, showering and transfers to the toilet.

The unique Hanger Bar features an Electrical Positioning System (EPS) and Anti-Rotation Mechanism (ARM) that allows caregivers to be in direct contact with the resident, providing an additional feeling of safety. The EPS also contributes to easy repositioning from sitting into a lying position, simplifying transfers making them effortless and more accurate. Within the rotatable hanger bar, the resident is given ample space to feel comfortable and secure without touching legs or knees, all while being lifted into the desired care position.

CARLO Comfort is available with a safe working load of 185 kg (408 lb) and 230 kg (507 lb) for bariatric residents with low mobility. Also available with low legs to accommodate lower beds, stretchers and furniture.

Battery and charger are backwards compatible across all Handicare by Beka mobile lifts.

# Caregiver and Resident Benefits



Carrybar with anti-rotation mechanism (ARM) and a free rotation area of 360°.



Electrically spreadable legs with low height to pass underneath baths and beds.



Four high-quality double castors, two of them with brakes, provide easy movement.



Electrical positioning system (EPS) in order to position your patients without any effort and force.



Cable remote control with display of the charging condition of the battery, failure diagnostics, service.



Light construction. Chassis, lifting column and components made of 100% aluminum.



Electrical height adjustment with lowering of the carrybar to the floor. Also includes manual emergency lowering system.



Wide choice of 4-point slings with safety clip system.

# CARLO Comfort Floor Lift

## Specifications

### Safe Working Load (SWL)

CARLO Comfort: 185 kg (408 lb)  
CARLO Comfort 230: 230 kg (507 lb)

### Dimensions (L x W x H)

#### CARLO Comfort

Carrier frame, length: 1290 mm (50.8 in)  
Carrier frame, height: 125 mm (4.9 in)

Inner width, closed: 540 mm (21.3 in)  
Inner width, spread: 1040 mm (40.9 in)

Outer width, closed: 700 mm (27.6 in)  
Outer width, spread: 1180 mm (46.5 in)

Total width of lift, carrier frame closed:  
700 mm (27.6 in)

Height: 1340 – 2170 mm (52.8 – 85.4 in)

#### CARLO Comfort 230

Carrier frame, length: 1400 mm (55.1 in)  
Carrier frame, height: 125 mm (4.9 in)

Inner width, closed: 540 mm (21.3 in)  
Inner width, spread: 1090 mm (42.9 in)

Outer width, closed: 700 mm (27.6 in)  
Outer width, spread: 1230 mm (48.4 in)

Total width of lift, carrier frame closed:  
700 mm (27.6 in)

Height: 1340 – 2170 mm (52.8 – 85.4 in)

For additional technical data, please refer to the user manual.

### Lifting Range

1400 mm (55.1 in)

### Castors

100 mm (3.9 in)

### Power Supply

Removable 24V battery

### Product Weight Without Packaging

55 – 65 kg (121 – 143 lb)

## Part Numbers

**925006008:** CARLO Comfort  
On board charging and dual controls

**925006108:** CARLO Comfort with Low Legs  
On board charging and dual controls

**925006018:** CARLO Comfort  
Includes two batteries, wall mount charger  
and dual controls

**925007008:** CARLO Comfort 230  
On board charging and dual controls

**925007018:** CARLO Comfort 230  
Includes two batteries, wall mount charger  
and dual controls

## Accessories / Options

### CARLO Comfort

Description	Part No.
Conversion kit for retrofitting a patient scale	924003400
Calibrated patient scale	924003300

### CARLO Comfort 230

Description	Part No.
Conversion kit for retrofitting a patient scale	924001800
Calibrated patient scale up to 250 kg (551 lb)	925001208

### For all floor lifts

Description	Part No.
Spare battery, 24V, blue	273225
Spare battery, 24V, orange	921070602

Are you ready to make safe  
and effortless resident transfers  
part of your everyday?

To learn more, call 1-888-637-8155 or  
visit [www.handicare.ca](http://www.handicare.ca) / [handicareusa.com](http://handicareusa.com)

### HANDICARE

485 Millway Avenue, Unit 2, Concord, ON, L4K 3V4

Toll Free: 877.304.5438 | Fax: 416.260.0256

[customerservice.canada@handicare.com](mailto:customerservice.canada@handicare.com)

[www.handicare.ca](http://www.handicare.ca)

10888 Metro Court, St. Louis, MO 63043

Toll Free: 866.891.6502 | Fax: 314.692.7858

[customerservice.northamerica@handicare.com](mailto:customerservice.northamerica@handicare.com)

[www.handicareusa.com](http://www.handicareusa.com)