

## Carina350 Dimensions

Name	Carina350EM		Carina350EE	
Part Number	60600011	60600012	60600009	60600014
Description	Manual base widening	Manual base widening with low legs	Electrical base widening	Electrical base widening with low legs
Material	Powder-coated steel			
Unit	1 piece	1 piece	1 piece	1 piece
Total height	50.8in – 72.4in (1290-1840mm)			
Base length	45.3in (1150mm)			
Base width (outer)	26.8in – 42.5in (680-1080mm)			
Base width (inner)	22in – 37in 560-940mm)			
Sling bar width	17.7in (450mm)			
Lifting speed	0.9in/sec			
Lifting range	22.8in – 66.1in (580-1680mm)			
Castors	3.9in (100mm) (front & rear)	1.4 in (front) 3.9 in (rear)	3.9in (100mm) (front & rear)	1.4 in (front) 3.9 in (rear)
Safe Working Load	350lb (159kg)			
Overall weight of lift	76lb (35kg)	78lb (35.4kg)	76lb (35kg)	78lb (35.4kg)