

## Eva450/Eva600 Dimensions

Name	Eva450EE	Eva450EE, low legs	Eva600EE	Eva600EE, low legs
Part Number	60100002	60100006	60100003	60100010
Description	Electrical base widening	Low legs with electrical base widening	Electrical base widening	Low legs with electrical base widening
Unit	1 piece	1 piece	1 piece	1 piece
Base width, outer dimensions	27in – 38.6in (685-80mm)		29.5in – 45.2in (750-1150mm)	
Base width, inner dimensions	23in – 34.6in (585-880mm)		25.6in – 41.3in (650-1050mm)	
Carry bar width	17.7in (450mm)		23.6in (600mm)	
Lifting speed	1.4in/sec (35mm/sec)			
Lifting range	19.7in – 70.5in (500-1790mm)		19.7in – 71.2in (500-1810mm)	
Castors	3.9in (100mm) front	1.4in (36mm) front	3.9in (100mm) front	1.4in (36mm) front
	3.9in (100mm) rear	3.9in (100mm) rear	3.9in (100mm) rear	3.9in (100mm) rear
Safe Working Load	450lb (205kg)		600lb (270kg)	
Overall weight of lift (including batteries)	85.9lb (39kg)		85.9lb (39kg)	