



Manual Transfer Aids

Positioning

Sliding Aids



EasyGlide 5020

EasyGlide

For sitting transfers from wheelchairs to beds, toilets, car seats and many other seated positions. It is especially suitable for leveling out height differences or for bridging short distances. Available in two models and several sizes.

Part No.	Product Name	Description	Unit
5020	EasyGlide	23.6 x 13.0 in (600 x 330 mm), three-part	1 Pc
5050	EasyGlide	29.5 x 13.0 in (750 x 330 mm), three-part	1 Pc
5060	EasyGlide	23.6 x 9.8 in (600 x 250 mm)	1 Pc
5070	EasyGlide	17.7 x 7.9 in (450 x 200 mm)	1 Pc



EasyGlide 5060

Accessories

Part No.	Product Name	Description	Unit
4014	Storage pocket	For EasyGlides	1 Pc



EasyGlide 5000

EasyGlide

For supine transfers from one bed to another.

Part No.	Product Name	Description	Unit
5000	EasyGlide, extra long	47.2 x 23.6 in (1200 x 600 mm)	1 Pc



EasyGlide 5045

EasyGlide

EasyGlide are used in pairs to facilitate positioning further back in chair. The smaller size makes it easier to place the sliding boards under the user's thighs and seat because they split the weight over the surface. EasyGlide, oval can also be used to facilitate the application of leg supports of slings with divided leg supports in a seated position. The sliding boards are also useful for helping a user lie down in bed. EasyGlides are easy to clean and wipeable.

Part No.	Product Name	Description	Unit
5045	EasyGlide, oval	17.7 x 7.9 in (450 x 200 mm)	2 Pcs





EasyGlide 5042

EasyGlide

For supine transfers from one bed to another. It is also for positioning higher up in bed and placing a sling behind the backrest of the user, wheelchair, chair or on the floor.

Part No.	Product Name	Description	Unit
5040	EasyGlide, oval	33.5 x 15.4 in (850 x 390 mm)	1 Pc
5042	EasyGlide, oval	33.5 x 15.4 in (850 x 390 mm)	2 Pcs



OneWaySlide 2075

OneWaySlide

OneWaySlide has many functions. The main one improves sitting positions in chairs or wheelchairs. The sliding mat facilitates movement in a backward direction, but will not slide forward, which makes it easier for the user to get and stay into a more comfortable position in a chair. OneWaySlide is also available in a tube-shaped design without handles for a more active user.



OneWaySlide 2069

Part No.	Product Name	Description	Unit
2074	OneWaySlide	43.3 x 11.8 in (1100 x 300 mm), 6 handles	1 Pc
2075	OneWaySlide	51.2 x 15.7 in (1300 x 400 mm), 6 handles	1 Pc
2076	OneWaySlide	51.2 x 17.7 in (1300 x 450 mm), 6 handles	1 Pc
2077	OneWaySlide	51.2 x 21.7 in (1300 x 550 mm), 6 handles	1 Pc
2069	OneWaySlide, tubular	19.7 x 15.7 in (500 x 400 mm)	1 Pc
2070	OneWaySlide, tubular	19.7 x 17.7 in (500 x 450 mm)	1 Pc



Positioning

Sliding Aids



EasySlide 1220

EasySlide

EasySlide (1000 x 700 mm) is mainly designed for users who need assistance to turn and move in bed. It performs best on a rough surface. EasySlide reduces the friction on both shoulder and hip, which reduces the risk of pressure sores. These features make EasySlide very gentle for users in pain.

Part No.	Product Name	Description	Unit
1220	EasySlide	39.4 x 27.6 in (1000 x 700 mm)	1 Pc
1230	EasySlide (1220 + 1251)	39.4 x 27.6 in (1000 x 700 mm), complete	1 Set



Disposable cover

Accessories

Part No.	Product Name	Description	Unit
1251	Protective cover	For 1220, 1230	1 Pc
1255	Disposable cover	For 1220, 1230	2 x 90 Pcs
1085	WallRack	For storing of sliding mats and disposable covers for EasySlide	1 Pc



WallRack



EasySlide

EasySlide (1400 x 900 mm) is designed for users with limited to severely reduced abilities to participate in the turning process. These users are often permanently confined to the bed, and using EasySlide helps to prevent pressure sores.

EasySlide (2000 x 900 mm) has extra padding in order to further reduce the risks of pressure sores. The padding distributes the pressure over the entire body and reduces pressure, friction and the risk of shearing on the especially vulnerable pressure points: the shoulders, hips, head and heel.



EasySlide 1430

Part No.	Product Name	Description	Unit
1420	EasySlide	55.1 x 35.4 in (1400 x 900 mm), 2 fastening straps	1 Pc
1430	EasySlide (1420 + 1451)	55.1 x 35.4 in (1400 x 900 mm), complete	1 Set
1470	EasySlide	78.7 x 35.4 in (2000 x 900 mm), 4 fastening straps	1 Pc
1480	EasySlide (1470 + 1481)	78.7 x 35.4 in (2000 x 900 mm), complete	1 Set



EasySlide 1470

Accessories

Part No.	Product Name	Description	Unit
1451	Protective cover	For 1420, 1430	1 Pc
1481	Protective cover	For 1470, 1480	1 Pc
1599	LockingStrap	Fastening straps, spare part for 1420, 1430, 1470, 1480	1 Pair
1085	WallRack	For storing of sliding mats and disposable covers for EasySlide	1 Pc




WallRack





EasySlide 2010



EasySlide 2015



Disposable cover



WallRack

EasySlide

EasySlide (500 x 600 mm), (500 x 700 mm) and (250 x 600 mm) are excellent aids for turning in bed, for moving users higher up in bed or for transferring them from the bed to the wheelchair. EasySlide is a small, practical aid that encourages activation. For the caregiver, EasySlide decreases the number of assisted positionings and eliminates unnecessary lifting, which spares their back and shoulders.

Part No.	Product Name	Description	Unit
2010	EasySlide	19.7 x 23.6 in (500 x 600 mm), polyester	1 Pc
2020	EasySlide	19.7 x 23.6 in (500 x 600 mm), microfiber	1 Pc
2030	EasySlide (2020 + 2051)	19.7 x 23.6 in (500 x 600 mm), complete	1 Set
2021	EasySlide	19.7 x 27.6 in (500 x 700 mm), microfiber	1 Pc
2015	EasySlide	9.8 x 23.6 in (250 x 600 mm), polyester	1 Pc

Accessories

Part No.	Product Name	Description	Unit
2051	Protective cover	For 2010, 2020, 2030	1 Pc
2052	Protective cover	For 2011, 2021	1 Pc
2055	Disposable cover	For 2010, 2020, 2030	100 Pcs
1085	WallRack	For storing of sliding mats and disposable covers for EasySlide	1 Pc



EasySlide, Disposable 2040



EasySlide, Disposable

EasySlide, disposable is a thin, pliable and very easy-to-use personal sliding mat for short-term use. EasySlide, disposable can be used to facilitate turning and repositioning higher up in bed, for supine and seated transfers, and for repositioning in the seated position. EasySlide, disposable can also be used to lift the user's legs in and out of bed, to assist dressing and undressing, and to facilitate application of a lifting sling.

EasySlide, disposable is a new, unbelievably versatile sliding mat that makes many common transfers, even heavy ones, easier for both the user and the caregiver. The sliding mat is made of environmentally friendly polyethylene (PE) plastic and the practical, hygienic package contains a roll with up to 60 individual EasySlide mats for short-term use.

Part No.	Product Name	Description	Unit
2040	EasySlide, disposable	19.7 x 27.6 in (500 x 700 mm)	2 x 60 Pcs





WendyLett

A comfortable functional material

The WendyLett sheets are made of a mixed material fabric that includes cotton, polyester and satin. Smooth cotton provides for comfort. Strong polyester ensures functionality and durability. Satin enables turning with minimal effort. The fabric quality is continuously monitored and evaluated for its strength, durability, burls, washing properties and non-inflammability.

Communicating and self-instructing design

The colors and patterns of the fabric in the WendyLett sliding and drawsheets are designed to indicate the function of the sheet and thereby also the needs of the user. White and grey stripes and checked patterns communicate which direction sliding is facilitated and how the sheet should be used. The design also makes the sheets easy to distinguish in the laundry. WendyLett sheets are also aesthetically pleasing.

A transfer aid which prevents pressure-related sores

By using WendyLett, small transfers, turning and positioning of users in bed can be performed without friction and shearing damage of the user's skin. The sheets are made of a fine material and can favorably be placed on top of a pressure-relieving mattress without compromising the pressure-relieving properties of the mattress. In order to provide the user with maximum comfort, all seams have been eliminated across the sheets and the joint seams are made of different functional materials that are flat and low.

Complete WendyLett Packages

Part No.	WendyLett, Sliding Sheet	WendyLett2Way, Sliding Sheet	Combination
1615	55.1 x 78.7 in (1400 x 2000 mm), white	35.4 x 70.9 in (900 x 1800 mm), white	(1614 + 1620)
1630	39.4 x 78.7 in (1000 x 2000 mm), grey / white striped	39.4 x 78.7 in (1000 x 2000 mm), grey / white striped	(1629 + 1631)
1635	55.1 x 78.7 in (1400 x 2000 mm), grey / white striped	55.1 x 78.7 in (1400 x 2000 mm), grey / white striped	(1634 + 1636)
Part No.	WendyLett, Sliding Sheet	WendyLett2Way, Sliding Sheet	Combination
1638	55.1 x 78.7 in (1400 x 2000 mm), grey / white striped	55.1 x 78.7 in (1400 x 2000 mm), grey / white checked	(1634 + 1637)
1628	55.1 x 78.7 in (1400 x 2000 mm), grey / white striped	55.1 x 74.8 in (1400 x 1900 mm), grey / white checked, 2 x 7 handles	(1634 + 1657*)
1644	78.7 x 35.4 in (2000 x 900 mm), grey / white striped, tailored with drawstring	55.1 x 78.7 in (1400 x 2000 mm), grey / white checked	(1641 + 1637)
1645	78.7 x 35.4 in (2000 x 900 mm), grey / white striped, tailored with drawstring	78.7 x 78.7 in (2000 x 2000 mm), grey / white checked	(1641 + 1647)

*Can be used in combination with mobile or stationary lifts.



WendyLett

WendyLett is intended for users who have difficulties turning in bed. WendyLett has a smooth, silky central part, enabling turning with minimal effort. The WendyLett base sheet is a satin sliding sheet. The sliding takes place sideways along the direction of the stripes. Sliding sideways is made easy without risk of sliding off the bed by the anti-slip material at both sides of the sheet. The material also prevents the user from sliding down in bed. WendyLett is used on the bed like a normal sheet with the smooth surface in the middle of the bed. For users who need more assistance, WendyLett can be combined with WendyLett2Way drawsheet or WendyLett4Way sliding drawsheet.



WendyLett, Sliding Sheet, White

Part No.	Description	Unit
1609	39.4 x 78.7 in (1000 x 2000 mm), sliding surface 23.6 in (600 mm), fits 31.5–35.4 in (800–900 mm) bed	1 Pc
1614	55.1 x 78.7 in (1400 x 2000 mm), sliding surface 23.6 in (600 mm), fits 31.5–35.4 in (800–900 mm) bed	1 Pc

WendyLett, Sliding Sheet, Grey / White Striped

Part No.	Description	Unit
1629	39.4 x 78.7 in (1000 x 2000 mm), sliding surface 23.6 in (600 mm), fits 31.5–35.4 in (800–900 mm) bed	1 Pc
1634	55.1 x 78.7 in (1400 x 2000 mm), sliding surface 23.6 in (600 mm), fits 31.5–35.4 in (800–900 mm) bed	1 Pc
1639	78.7 x 78.7 in (2000 x 2000 mm), sliding surface 27.6 in (700 mm), fits 41.3 in (1050 mm) bed	1 Pc
1667	39.4 x 86.6 in (1000 x 2200 mm), sliding surface 39.4 in (1000 mm), fits 47.2 in (1200 mm) bed	1 Pc
1668	55.1 x 86.6 in (1400 x 2200 mm), sliding surface 39.4 in (1000 mm), fits 47.2 in (1200 mm) bed	1 Pc
1669	78.7 x 86.6 in (2000 x 2200 mm), sliding surface 39.4 in (1000 mm), fits 47.2 in (1200 mm) bed	1 Pc
1681	78.7 x 143.7 in (2000 x 3650 mm), sliding surface 27.6 in (700 mm), extra wide to allow for optional placing of sliding surface	1 Pc



WendyLett, Sliding Sheet, Grey / White Striped, Fitted Sheet with Elastic

Part No.	Description	Unit
1641	78.7 x 35.4 in (2000 x 900 mm), sliding surface 27.6 in (700 mm), fits 36 in (900 mm) bed	1 Pc
1651	78.7 x 41.3 in (2000 x 1050 mm), sliding surface 31.5 in (800 mm), fits 42 in (1050 mm) bed	1 Pc
1661	78.7 x 47.2 in (2000 x 1200 mm), sliding surface 37.4 in (950 mm), fits 47.2 in (1200 mm) bed	1 Pc
1671	78.7 x 55.1 in (2000 x 1400 mm), sliding surface 43.3 in (1100 mm), fits 55.1 in (1400 mm) bed	1 Pc



Positioning


Sliding Aids



WendyLett2Way

WendyLett2Way is a drawsheet made of a material which is adapted to WendyLett base sheet so that sliding down in bed is significantly reduced. Sliding takes place sideways, along the direction of the stripes. WendyLett2Way is placed on top of a WendyLett base sheet. The combination makes it easy for the caregiver to turn the user sideways.

Part No.	Product Name	Description	Unit
1620	WendyLett2Way	Drawsheet, 35.4 x 70.9 in (900 x 1800 mm), white	1 Pc
1631	WendyLett2Way	Drawsheet, 39.4 x 78.7 in (1000 x 2000 mm), white / grey striped	1 Pc
1636	WendyLett2Way	Drawsheet, 55.1 x 78.7 in (1400 x 2000 mm), white / grey striped	1 Pc



WendyLett4Way

WendyLett4Way is a sliding drawsheet to be combined with WendyLett base sheet. It offers multidirectional sliding, as indicated by the checked pattern. WendyLett4Way is placed on top of a WendyLett base sheet. The combination allows for multidirectional moves and facilitates procedures like turning in bed, moving higher up/lower in bed or other re-adjusting of a lying user, even with heavier or passive users.

Part No.	Product Name	Description	Unit
1657*	WendyLett4Way, with handles	Sliding drawsheet with 7 handles on each side, 55.1 x 74.8 in (1400 x 1900 mm), grey / white checked	1 Pc
1637	WendyLett4Way	Sliding drawsheet, 55.1 x 78.7 in (1400 x 2000 mm), grey / white checked	1 Pc
1647	WendyLett4Way	Sliding drawsheet, 78.7 x 78.7 in (2000 x 2000 mm), grey / white checked	1 Pc
1648	WendyLett4Way, inco	Sliding drawsheet / incontinence protection with extra fabric on both sides, 29.5 x 78.7 in (750 x 2000 mm), grey / white checked	1 Pc
1649	WendyLett4Way, inco with handles	Sliding drawsheet / incontinence protection with 7 handles on each side, 55.1 x 47.2 in (1400 x 1200 mm), grey	1 Pc
1652	WendyLett4Way	Sliding drawsheet, 55.1 x 102.4 in (1400 x 2600 mm), grey / white checked, for heavier users	1 Pc

*Can be used in combination with mobile or stationary lifts





WendyLean

The patented, self-retracting handles are handy when you need to help a user sit up in bed. The underside of the pillowcase is made of smooth satin which makes it easy to pull the pillow in under the shoulders for raising the upper body of the user.

Part No.	Product Name	Description	Unit
1680	WendyLean	19.7 x 23.6 in (500 x 600 mm), grey / white striped	1 Pc



Patent: SE-0400590-6



Positioning

Turning Plates



TurnTableSoft 5086

TurnTableSoft

TurnTableSoft makes it easier to turn when sitting down. TurnTableSoft is constructed of two round discs with a friction-reducing material on the inside, enabling the discs to rotate smoothly over each other.

Part No.	Product Name	Description	Unit
5086	TurnTableSoft	17.7 in dia (Ø 450 mm)	1 Pc
5087	TurnTableSoft	15.7 in dia (Ø 400 mm)	1 Pc



Other Positioning Aids



SitWalk 6080

SitWalk

This device is named after the method of moving a user back-wards in the chair to a more comfortable sitting position. Lean the user sideways and move one buttock at a time, facilitating a hip-hitch.

Part No.	Product Name	Description	Unit
6080	SitWalk	55.1 x 5.1 in (1400 x 130 mm), 8 handles	1 Pc



AntiSlip

AntiSlip comes in handy when you need a good grip for moving higher up in bed. It can also be used as a restraining surface between materials in many of our sliding products.

Part No.	Product Name	Description	Unit
8040	AntiSlip	17.7 x 15 in (450 x 380 mm)	1 Pc
8045	AntiSlip	17.7 x 15 in (450 x 380 mm)	10 Pcs
8046	AntiSlip	17.7 x 15 in (450 x 380 mm)	30 Pcs



FootStool 8090



FootStool 8095

FootStool

FootStool is used to raise the position of the user's knees to enable an easier and safer transfer in a seated position. FootStool also provides excellent support under the user's feet during visits to the toilet, or when repositioning further back in a chair or wheelchair.

Part No.	Product Name	Description	Unit
8090	FootStool	12.6 x 15 x 3.9 in (320 x 380 x 100 mm)	1 Pc
8095	FootStool	12.6 x 15 x 2 in (320 x 380 x 50 mm)	1 Pc



Design protection: RCD-67533-1



FootStoolCover

FootStoolCover is a wipeable cover designed to be used with the FootStool. It is made of wipeable material, which provides for optimal hygiene. The cover has a handle to facilitate transportation of the FootStool, or it can be used to hang up the FootStool when not in use.



Part No.	Product Name	Description	Unit
8091	FootStoolCover	Wipeable cover to FootStool (Part No. 8090)	1 Pc
8096	FootStoolCover	Wipeable cover to FootStool (Part No. 8095)	1 Pc



Design protection: RCD-67533-1

Support

Belts



FlexiBelt



FlexiBelt, hug



EasyBelt



EasyBelt, hug

FlexiBelt / EasyBelt

FlexiBelt and EasyBelt are used for sit-to-stand transfers or moving from the wheelchair to the bed, toilet or car. FlexiBelt has a soft inner side of polyester velvet that makes it comfortable to wear. EasyBelt has an anti-slip inside of PU-coated nylon that can be used in environments where moisture and dampness are present or when there are special hygiene requirements. EasyBelt is also available in a wipeable material that can be cleaned using disinfectants. FlexiBelt, hug and EasyBelt, hug are new versions of Handicare's classic belts. FlexiBelt, hug and EasyBelt, hug feature a new patented* design with horizontal handles that slide freely in loops to make the belt "hug" onto the user instead of sliding up during transfer situations. With FlexiBelt, hug and EasyBelt, hug the caregiver can provide even more support without using more effort. The new design does all the work by itself.

FlexiBelt

Part No.	Description	Handles	Unit
6023	S: 23.6–39.4 in (600–1000 mm)	7	1 Pc
6024	M: 27.6–47.2 in (700–1200 mm)	9	1 Pc
6025	L: 39.4–63.0 in (1000–1600 mm)	11	1 Pc
6026	XL: 45.3–68.9 in (1150–1750 mm)	13	1 Pc

FlexiBelt, hug

Part No.	Description	Handles	Unit
6072	XS: 17.7–27.6 in (450–700 mm)	3 + 2 movable	1 Pc
6073	S: 23.6–39.4 in (600–1000 mm)	4 + 3 movable	1 Pc
6074	M: 27.6–47.2 in (700–1200 mm)	5 + 4 movable	1 Pc
6075	L: 39.4–63 in (1000–1600 mm)	6 + 5 movable	1 Pc
6076	XL: 45.3–68.9 in (1150–1750 mm)	7 + 6 movable	1 Pc

EasyBelt

Part No.	Description	Handles	Unit
6012	XS: 17.7–27.6 in (450–700 mm)	5	1 Pc
6013	S: 23.6–39.4 in (600–1000 mm)	7	1 Pc
6014	M: 27.6–47.2 in (700–1200 mm)	9	1 Pc
6015	L: 39.4–63.0 in (1000–1600 mm)	11	1 Pc

EasyBelt, hug

Part No.	Description	Handles	Unit
6002	XS: 17.7–27.6 in (450–700 mm)	3 + 2 movable	1 Pc
6003	S: 23.6–39.4 in (600–1000 mm)	4 + 3 movable	1 Pc
6004	M: 27.6–47.2 in (700–1200 mm)	5 + 4 movable	1 Pc
6005	L: 39.4–63.0 in (1000–1600 mm)	6 + 5 movable	1 Pc
6006	XL: 45.3–68.9 in (1150–1750 mm)	7 + 6 movable	1 Pc

Accessories

SupportStraps

Part No.	Description	Unit
6061	Accessory for belts	1 Pair

CE *Patent: SE-536273, EP-2611409, Patent pending: JP-2013-527040



ReTurnBelt

ReTurnBelt

Our patented* ReTurnBelt is an excellent complement to ReTurn when there is a need for greater support and stability, including when the user's weight-bearing ability is questionable or varies. ReTurnBelt can also be used separately, for support during sit-to-stand and to stabilize the user in an upright position, such as during standing training, gait training and transfer. ReTurnBelt has an anti-slip inside of PU-coated nylon that can be used in environments where moisture and dampness are present or when there are special demands on hygiene. The fabric also ensures that the belt stays in place. ReTurnBelt is also available in a wipeable material that can be cleaned using disinfectants as well as a disposable material that can help prevent the spread of multi-resistant and other types of contagious bacteria.

ReTurnBelt

Part No.	Description	Handles	Unit
6031	XXS: 25.6–31.5 in (650–800 mm)	4	1 Pc
6032	XS: 31.5–35.4 in (800–900 mm)	4	1 Pc
6033	S: 35.4–41.3 in (900–1050 mm)	4	1 Pc
6034	M: 41.3–49.2 in (1050–1250 mm)	4	1 Pc
6035	L: 49.2–57.1 in (1250–1450 mm)	4	1 Pc
6036	XL: 53.2–61.0 in (1350–1550 mm)	4	1 Pc
6037	XXL: 61.0–68.9 in (1550–1750 mm)	4	1 Pc
6038	XXXL: 68.9–76.8 in (1750–1950 mm)	4	1 Pc

ReTurnBelt,
disposable

ReTurnBelt, disposable, non-woven

Part No.	Description	Handles	Unit
6233	S: 35.4–41.3 in (900–1050 mm)	2	5 Pcs
6234	M: 41.3–49.2 in (1050–1250 mm)	2	5 Pcs
6235	L: 49.2–57.1 in (1250–1450 mm)	2	5 Pcs
6236	XL: 53.2–61.0 in (1350–1550 mm)	2	5 Pcs



*Patent: SE-0401517-8

Support

Belts / Walking Aids



SupportStraps

SupportStraps is a multi-functional accessory providing extra support and security for users during standing, gait training and transfers with Handicare's belts, sit-to-stand aids and transfer platforms.

SupportStraps can be used as:

- Leg harness for EasyBelt / FlexiBelt
- Torso / side support when using QuickMove together with EasyBelt / FlexiBelt / FlexiMove
- Side support when using ReTurn together with ReTurnBelt

Part No.	Product Name	Description	Unit
6061	SupportStraps	Accessory for belts	1 Pair



DisposableBelt

DisposableBelt is intended for short-term use, as extra support during sit-to-stand transfers and sitting transfers from the wheelchair to the bed, toilet or car. DisposableBelt can also be used as support when the user is standing or walking.

Part No.	Product Name	Description	Unit
6054	DisposableBelt	One size, medium / large	5 Pcs



WalkingBelt

WalkingBelt is used to support children with impaired strength and stability during gait training. WalkingBelt, which provides support around the torso, is held in place with a padded leg harness. The length of the two padded handles can be adapted to give the caregiver a comfortable working position when walking close to the child holding the handles to provide the child with additional support and stability during gait training.

Part No.	Product Name	Description	Unit
6041	WalkingBelt	XXS: 23.6 in (600 mm)	1 Pc
6042	WalkingBelt	XS: 27.6 in (700 mm)	1 Pc



Find the Right Transfer Aid for the Right Situation

Handicare offer solutions for all transfer situations and needs. To make it easier to find and choose the right transfer aid for the right situation, please use the matrix below to combine the main products in the transfer categories with the 14 most common activities in user transfers.

PRODUCTS	ACTIVITIES													
	Supine Transfer	Turning in Bed	Higher Up in Bed	Lying to Sitting	Support for Chair Training	Standing Training	From Bed to Wheelchair	From Bed to Commode	From Wheelchair to Toilet	Raising – Sitting Down	Positioning in Chair	Into and Out of Vehicles	Positioning in Car	Preventing Pressure Sores
POSITIONING														
EasyGlide	✓						✓	✓	✓		✓	✓		
One Way Slide											✓			
EasySlide / EasySlide Disposable	✓	✓	✓				✓				✓			✓
WendyLett		✓												✓
WendyLett2Way		✓		✓										✓
WendyLett4Way		✓	✓	✓										✓
WendyLean	✓	✓	✓	✓										✓
Turn TableSoft				✓								✓		
SitWalk			✓								✓			
AntiSlip									✓					
FootStool / FootStoolCover							✓	✓	✓		✓			
SUPPORT														
FlexiGrip				✓			✓							
EasyBelt / FlexiBelt					✓		✓	✓	✓	✓	✓			
Disposable Belt / Walking Belt					✓		✓	✓	✓	✓	✓			



HANDICARE

10888 Metro Court, St. Louis, MO 63043

Toll Free: 866.891.6502 | Fax: 314.692.7858

customerservice.northamerica@handicare.com

www.handicareusa.com

485 Millway Avenue, Unit 2, Concord, ON, L4K 3V4

Toll Free: 877.304.5438 | Fax: 416.260.0256

customerservice.northamerica@handicare.com

www.handicare.ca

Learn more about safe
and easy transfer aids today!

Call 1-888-637-8155 or visit
www.handicareusa.com / handicare.ca