

ReTurn Dimensions

Name	ReTurn 7400	ReTurn7500i*	ReTurn7600
Part Number	7400	7500i	7600
Unit	1 piece	1 piece	1 piece
Baseplate Length	27in (685mm)	27in (685mm)	31.3in (795mm)
Baseplate Width	22.4in (570mm)	22.4in (570mm)	26.5in (672mm)
Baseplate Height	1.6in (40mm)	1.6in (40mm)	1.8in (45mm)
Overall Height	39.2in (995mm)	45.3in (1150mm)	47.6in (1210mm)
Safe Working Load	330lb (150kg)	330lb (150kg)	450lb (205kg)
Total weight of lift	36.4lb (16.5kg)	36.4lb (16.5kg)	57.3lb (26.0kg)

Patent: SE-0300528-7

*Patent pending: EP-12868688-8; US-14/378 658

Design Protection: RCD-529631-1