

Pressure Fit Lift System Dimensions

Name	1-Post				
Part Number	341525	341527			
Description	Wall mount with swivel trolley	Wall mount with standard trolley			
Unit	1 system	1 system			
Height	Expands from 7ft – 9ft				
Safe Working Load	440lbs				

Name	2-Post						
Part Number	341500	341503	341520	341521	341522	341526	341528
Description	2-Post	With C-300 trolley	Bath mounted	Bath mounted with swivel trolley	Bath mounted with C-300 trolley	Wall mounted with swivel trolley	Wall mounted with standard trolley
Unit	1 system	1 system	1 system	1 system	1 system	1 system	1 system
Adjustable width	7ft – 10ft						
Unit weight	45lbs						



Name	3-Post					
Part Number	341505	341506	341507	341515	341516	341517
Description	3-Post	With swivel trolley	With C-300 trolley	Bath mounted	Bath mounted with swivel trolley	Bath mounted with C-300 trolley
Unit	1 system	1 system	1 system	1 system	1 system	1 system
Adjustable width	6ft - 8ft					
Adjustable length	7ft – 10ft					
Unit weight	65lbs					

Name	4-Post				
Part Number	341510	341511	341512		
Description	4-Post	With swivel trolley	With C-300 trolley		
Unit	1 system	1 system	1 system		
Adjustable width		6ft – 8ft			
Adjustable length	7ft – 10ft				
Unit weight	100lbs				