



ANA Nursing and Care Table

Compact and Comfortable Universal Care



All-Purpose Nursing Table Designed for **Everyday Care**

ANA is designed to assist with patient care in a safe and comfortable manner and is ideal for changing places environments. The large lifting range allows for easy transfers from wheelchairs and optimum working heights for everyday care. The surface features moisture-proof, soft, comfortable padding that can easily be wiped down by caregivers for disinfection. The low profile of ANA makes it suitable for smaller rooms where space is a concern. Due to its wall mounted construction, no substructure exists that can harm children or adults when being used. A foldable model is also available allowing for additional floor space.


Caregiver and Resident Benefits



The included hand control allows for easily controlling the functions of ANA.



Additional accessories include side rails for roll-off protection and padding.



Moisture-proof, soft, comfortable padding for easy disinfection.



The electronic height adjustment and foldable version of ANA allows for ergonomic work in smaller, confined spaces such as bathrooms.



The stretcher can move down to 400 mm (15.7 in) for direct wheelchair transfers.



Choose from a wide selection of optional colored padding options.



ANA is available in different stretcher sizes to meet patient and room requirements. Optionally, it can be equipped with foldable side rails for extra safety.



Easily adjust ANA electronically for optimum working heights.



ANA is available in a foldable version making it perfect for rooms with limited space.

ANA Nursing and Care Table

Specifications

Part Numbers

Safe Working Load (SWL)

135 kg (298 lb)

Dimensions (L x W x H)

Length (optional):

1000 mm (39.4 in), 1400 mm (55.1 in),

1800 mm (70.9 in), 2000 mm (78.8 in)

Width (optional):

800 mm (31.5 in), 1000 mm (39.4 in)

Adjustable Height:

400–950 mm (15.7–37.4 in)

For additional technical data, please refer to the user manual.

Electrical Height Adjustment

Description	Part No.
ANA EL 100-80	950111000
ANA EL 140-80	950112000
ANA EL 180-80	950113000
ANA EL 200-80	950114000
ANA EL 140-100	950115000
ANA EL 180-100	950116000
ANA EL 200-100	950117000

Electrical Height Adjustment With Foldable Stretcher

Description	Part No.
ANA EL 100-80, foldable	950111001
ANA EL 140-80, foldable	950112001
ANA EL 180-80, foldable	950113001
ANA EL 200-80, foldable	950114001
ANA EL 140-100, foldable	950115001
ANA EL 180-100, foldable	950116001
ANA EL 200-100, foldable	950117001

Accessories / Options

Description	Part No.
Rail protection front bar	950117100
Padding for front bar	950117150
Roll protection head and foot side bars (set)	950117200
Padding for head and foot side bars (set)	950117250
Roll protection rear bars (set)	950117300
Padding for rear bars (set)	950117350
Head part	950109000

HANDICARE

81 Romina Drive, Concord, ON, L4K 4Z9

Toll Free: 877.304.5438 | Fax: 416.260.0256

customerservice.northamerica@handicare.com

www.handicare.ca

10888 Metro Court, St. Louis, MO 63043

Toll Free: 866.891.6502 | Fax: 314.692.7858

customerservice.northamerica@handicare.com

www.handicareusa.com

Are you ready to make daily care
safer and more comfortable?

To learn more, call 1-888-637-8155 or
visit www.handicare.ca / handicareusa.com