



MONA Care and Treatment Table

Designed for Comfort and Mobility



The Ideal Solution for **Everyday Care**

The MONA Care and Treatment Table is comprised of a large comfortable support pad designed for patients and residents who require assistance with everyday care. Designed with safety in mind, the table features rounded padded ends to prevent injury to exposed susceptible skin. The electronic height adjustment features a unique scissor mechanism that eliminates wobbling and reinforces stability for even the most restless patients, ensuring caregivers can work ergonomically and safely.

MONA can be furnished with a folding side rail on each side providing additional safety. A sturdy metal fastener ensures that the side rails cannot fold down, even under strong pressure. Continuous height adjustment (electric or hydraulic options) allows optimum work and transfer height to be set for care tasks. The base frame can be accessed underneath from all sides and is equipped with four wheels and locking castors for easy mobility.

Caregiver and Resident Benefits



Additional accessories include side rails for roll-off protection and padding.



Available in electronic or hydraulic height adjustment options.



Choose from a wide selection of optional colored padding options.



Easily integrates with floor lifts providing adequate clearance under the frame for access.



Scissor mechanism eliminates wobbling and reinforces stability.



Large comfortable support pad with rounded ends for increased safety.



MONA is available in different stretcher sizes to meet patient and room requirements. Optionally, it can be equipped with folding side rails for extra safety.



Easily adjust MONA electronically with the included hand control for optimum working heights.



MONA comes equipped with four wheels featuring locking castors for easy transport.

MONA Care and Treatment Table

Specifications

Safe Working Load (SWL)

Electric or hydraulic
205 kg (452 lb)

Dimensions (L x W x H)

Pad length (optional) – special sizes upon request:
1500–2000 mm (59.1–78.8 in)

Pad width (optional) – special sizes upon request:
650 mm (25.6 in), 800 mm (31.5 in), 1000 mm (39.4 in)

Lifting height (hydraulic height adjustment):
480–990 mm (18.9–39.0 in)

Lifting height (electrical height adjustment):
480–990 mm (18.9–39.0 in)

Clearance under frame
138 mm (5.4 in)

For additional technical data, please refer to the user manual.

Part Numbers

Electrical Height Adjustment

Description	Part No.
MONA EL 150-65	950100000
MONA EL 150-80	950102000
MONA EL 200-80	950104000
MONA EL 200-100	950106000

Hydraulic Height Adjustment

Description	Part No.
MONA HYD 150-65	950101000
MONA HYD 150-80	950103000
MONA HYD 200-80	950105000
MONA HYD 200-100	950107000

Accessories / Options

Description	Part No.
Rail protection side bar	950103400
Padding for side bar	950103410
Head part	950109000
Chassis cover made of ABS	950103300

HANDICARE

81 Romina Drive, Concord, ON, L4K 4Z9
Toll Free: 877.304.5438 | Fax: 416.260.0256
customerservice.northamerica@handicare.com
www.handicare.ca

10888 Metro Court, St. Louis, MO 63043
Toll Free: 866.891.6502 | Fax: 314.692.7858
customerservice.northamerica@handicare.com
www.handicareusa.com

Are you ready to make daily care
safer and more comfortable?

To learn more, call 1-888-637-8155 or
visit www.handicare.ca / handicareusa.com