



# Bure Rise & Go

Early Mobilization and Rehabilitation Aid



## Ergonomically designed, powered height adjustment, and frame widening walker with power rise

The Bure Rise & Go provides a patented power rise function combining stand-up, mobility, and a walking aid all in a single product. Rise & Go is equipped with powered frame widening and low leg height making it a versatile walker for wheelchair and bed applications in both acute hospital and long-term care environments.

# Caregiver and User Benefits



Ergonomic handles adjustable in any direction for optimum user comfort, added benefit for stroke patients who can reach their "usual" hand position.



Electrically driven height adjustment makes it easy to pre-set to the desired height. Once the user is standing in position the height can be adjusted for comfort and proper support for optimal walking.



All four castors are individually lockable as standard.



Hand control conveniently located for both user and caregiver to easily access. All Bure Rise & Go parts have IP65 international protection ratings and can therefore be used in wet spaces.



The power rise harness is designed to reach as far down across the bottom as possible, providing optimum support resulting in a more natural standing action.



Bure Rise & Go is fitted with shin supports that provide relief to the knees during sit-to-stand transfers, they are both height adjustable and can be moved aside to create greater space.



Anatomical pads provide relief for shoulders, arms, and neck allowing focus to be shifted to leg and chest muscles, improving stability, and making Bure easier to steer. The anatomically designed armrests can also be folded out to provide users with an easy means of reaching the handles during sit-to-stand transfers.



The open design provides stability with greater room to maneuver. Users do not feel shut in and caregivers are able to assist with standing up more easily. With the frame positioned outward, Bure becomes more stable and helps those who require extra walking space.



Electronic frame widening makes it easy to move Bure right up close to users without worrying about wheelchair width, chair width or bed design.

# Bure Rise & Go

## Specifications

### Safe Working Load (SWL)

150 kg (330.7 lb)

### Dimensions (Width)

660 mm – 106 mm (in)

### Dimensions (Length)

1010 mm (in)

### Dimensions (Height)

900 - 1300 mm (in)

### Dimensions (Armrest / Internal)

340 – 670 mm

## Part Numbers

**56-309/75:** Bure Rise & Go DB with electric frame widening, 75 mm

**56-309:** Bure Rise & Go DB with electric frame widening

## Accessories / Options

### Platform

Description	Part No.
Short	56-333-DB
Long	56-383-DB

### Stand-up Harness, Multi

Description	Part No.
Small	56-384-Multi-S
Medium	56-384-Multi-M
Large	56-384-Multi-L
X-Large	56-384-Multi-XL

### Gait Training Kit

Description	Part No.
For Multi Stand-up Harness	56-388-KIT

### HANDICARE

81 Romina Drive, Concord, ON, L4K 4Z9

Toll Free: 877.304.5438 | Fax: 416.260.0256

[customerservice.northamerica@handicare.com](mailto:customerservice.northamerica@handicare.com)

[www.handicare.ca](http://www.handicare.ca)

10888 Metro Court, St. Louis, MO 63043

Toll Free: 866.891.6502 | Fax: 314.692.7858

[customerservice.northamerica@handicare.com](mailto:customerservice.northamerica@handicare.com)

[www.handicareusa.com](http://www.handicareusa.com)

Are you ready to get started  
with the Bure Rise and Go?

To learn more, call 1-888-637-8155 or  
visit [www.handicare.ca](http://www.handicare.ca) / [handicareusa.com](http://handicareusa.com)