



# Encore with Retractable Deck Technology




**Encore:** seat deck retracts during HOB elevation maintaining the user in a safer comfortable zone of vertical alignment (shaded green).

Why vertical alignment is vital in long-term care:

- Superior fall prevention When the user is kept in vertical alignment all access-points in the user's bed environment remain safely within reach.
- Superior staff protection Fixed-deck beds cannot keep the user in alignment necessitating frequent turning or repositioning of the user, increasing the prospect of staff injury.



**Ordinary fixed deck bed:** the seat deck does not move during HOB elevation, forcing user out of vertical alignment.



## Advanced Clinical Positioning

The Encore foot-end control provides one-button advanced clinical positioning with staff-only lockout.

