

# Rexx and Fast Rexx

A favorite Q-Series design available in two distinct models



**Q-Series**  
Quality. Detail. Value.

- The Rexx which rises from 7-7/8" low to 26" high
- The fast-rising Fast Rexx which rises from 7-7/8" to 28" high in 20 seconds

The Fast Rexx meets OSHA speed recommendation of 20 seconds; can be expected to reduce staff injuries related to stooping and lifting, increase staff productivity and reduce staff turnover rates.


Standard Pendant Controller with Optional Foot-end Control



Item #	Description	Item #	Description
QD2000B	Rexx bed	Q6412LP	Rexx side rail mount pendant holder
QD2250B	Rexx Fast bed	QP6417LP	Rexx gooseneck pendant holder - side rail mount
QD2000M	Rexx bed with 360° mobility	Q6416LP	Rexx plastic pendant holder - head/foot board mount
QD2000ML	Rexx bed with mobility and floor locks	Q6417LP	Rexx gooseneck pendant holder - board mount
QD2000MLF	Rexx bed with mobility, floor locks and Foot End Controls	Q6730	Universal trapeze bar
QD2250ML	Rexx bed with mobility and floor locks	Q6733	Advantage trapeze mounting frame
QD2250MLF	Rexx bed with mobility, floor locks and Foot End Controls	Q6584	Fast Rexx pendant holder - side rail mount
RR938-L	Three position pivot assist - left	Q6591	Fast Rexx gooseneck pendant holder - side rail mount
RR938-R	Three position pivot assist - right	Q6589	Fast Rexx plastic pendant holder - head/foot board mount
RR939	Three position pivot assist - set	Q6590	Fast Rexx gooseneck pendant holder - board mount
RM930	Three position pivot assist rail - universal		
RM950	Fixed assist rail - universal		
RAD306H	Standard 1/2 head end side rail - set		
RR501	3/4 length collapsible side rails		

Height Floor to Mattress Support	7-7/8" - 26" (28" Fast Rexx)
Deck Angles - Back Elevation	0 - 70°
Deck Angles - Knee Elevation	0 - 20°
Bed Length - Overall Length	84-1/2" / 88-1/2"
Bed Length - Sleep Surface	76" / 80"
Bed Width - Sleep Surface	35"
Weight of Bed (without accessories)	170 lbs / 78 kg
Safe Working Load	450 lbs / 204 kg

Bed frame, overhead view



Convenient decals indicate proper placement of each expander.

## Tool-less 42" Expander Kits for Rexx and Advantage Beds

Easy on, easy off expander kits install instantly.

Require no tools, bolts, or pins.

Each kit contains 6 deck expanders (3 each right and left) required to expand a 35" bed to 42" wide.

Item #	Description
Q6709	42" Deck Expander Kit for Rexx, Fast Rexx, and Advantage models



Kit is supplied with protective packaging for storage and re-use.



For ordering and technical support on beds, accessories and furniture, call: **800-268-8671**