

Span F450T

Standing Transfer Aid

- For users requiring Minimum Assistance
- 450 lb. (205 kg) safe working load
- Keeps users actively engaged in the transfer process
- **Easy base width adjustment** accommodates chairs and commodes.
- **A convenient toileting alternative** vs. use of a wheelchair.
- **Locking rear casters** add an extra measure of safety.
- **Attachment hooks** for optional support belt.
- **Simple, fast assembly** using tools provided.



SPECIFICATIONS

| | |
|-----------------------------|-----------------|
| Weight Capacity | 450 lb / 205 kg |
| Overall Height | 43.2" |
| Overall Base Length | 35.5" |
| Minimum Seat Height | 26" |
| Maximum Seat Height | 33.5" |
| Front-end Base Height | 3.6" |
| Internal Leg Width (Open) | 32" |
| External Leg Width (Closed) | 27.7" |
| Horizontal Handlebar Width | 20" |
| Net Weight | 78 lb / 35.2 kg |
| Shipping Weight | 93 lb / 42 kg |

Item # Description

| | |
|----------|---------------------------|
| F450T | Sit-To-Stand Patient Lift |
| Warranty | 3 years (frame) |

Span F500P

Full Body Patient Lift

- For users requiring Maximum Assistance
- 500 lb. (228 kg) safe working load
- 6-point spreader bar for maximum sling options
- **Easy base width adjustment** with locked position; accommodates chairs and commodes.
- **Redundant operation** via control unit, allowing use if pendant is lost or damaged.
- **Floor extraction:** Can be used to lift patients from floor.
- **Overload protection** with visual/audible alert, halts lifting when capacity is exceeded.
- **Manual emergency lowering** for added safety.
 - **Quick-release battery** can be charged via charging cradle (optional), or via plug-in charger (included) for charging in place.
 - **Backlit LCD display** during operation for ease of viewing.
 - **Simple, fast assembly** using tools provided.



SPECIFICATIONS

| | |
|----------------------------|-----------------|
| Weight Capacity | 500 lb / 228 kg |
| Base Legs Length | 44.9" |
| Base Closed External Width | 24.9" |
| Base Open External Width | 42.9" |
| Mast Height | 48.3" |
| Base Height | 4.5" |
| Minimum Lifting Height | 25.8" |
| Maximum Lifting Height | 69.1" |
| Net Weight | 99 lb / 45 kg |
| Shipping Weight | 108 lb / 49 kg |

Item # Description

| | |
|----------|--|
| F500P | Full Body Patient Lift |
| Warranty | 3 years (frame and spreader bar); 2 years (electronics); 1 year (battery) |

