

# Span F500S

## Powered Sit-to-Stand Patient Lift

- For users requiring Moderate Assistance
- 500 lb. (228 kg) safe working load
- Contoured knee pad with three height settings
- **Easy base width adjustment** with locked position; accommodates chairs and commodes.
- **Wide non-slip footplate** provides stable base for transfer.
- **Redundant operation** via control unit, allowing use if pendant is lost or damaged.
- **Overload protection** with visual/audible alert, halts lifting when capacity is exceeded.
- **Manual emergency lowering** for added safety.
- **Quick-release battery** can be charged via charging cradle (optional), or via plug-in charger (included) for charging in place.
- **Backlit LCD display** during operation for ease of viewing.
- **Simple, fast assembly** using tools provided.



### SPECIFICATIONS

Weight Capacity	500 lb / 228 kg
External Base Width (Closed)	24.9"
Internal Leg Width (Open)	33.9"
Upright Mast Height	40"
Overall Base Length	37.5"
Base Height	4.9"
Minimum Lifting Height with Spreader Bar	38"
Maximum Lifting Height with Spreader Bar	60.5"
Knee Pad Length Extension Range	4.5"
Maximum Boom Width	17.7"
Net Weight	111 lb / 50.3 kg
Shipping Weight	119 lb / 53.9 kg

Item #	Description
F500S	Sit-To-Stand Patient Lift
Warranty	3 years (frame and spreader bar); 2 years (electronics); 1 year (battery)

# Span F600B

## Bariatric Full Patient Lift

- For users requiring Maximum Assistance
- 600 lb. (273 kg) safe working load
- 6-point spreader bar for maximum sling options
- **Locking rear casters**, foot-activated lock/unlock.
- **Floor extraction:** Can be used to lift patients from floor.
- **Quick-release battery** can be charged via charging cradle (optional), or via plug-in charger (included) for charging in place.
- **Overload protection** with visual/audible alert, halts lifting when capacity is exceeded.
  - **Manual emergency lowering** for added safety.
  - **Simple, fast assembly** using tools provided.



### SPECIFICATIONS

Weight Capacity	600 lb / 273 kg
Overall Base Length	46.9"
External Base Width (Closed)	26.4"
Internal Leg Width (Open)	39.3"
Upright Mast Height	51.3"
Base Height	4.7"
Minimum Lifting Height with Spreader Bar	21.5"
Maximum Lifting Height with Spreader Bar	68.6"
Net Weight	123 lb / 56 kg
Shipping Weight	139 lb / 63 kg

Item #	Description
F600B	Bariatric Full Body Patient Lift
Warranty	3 years (frame and spreader bar); 2 years (electronics); 1 year (battery)

