

# Span F400

## Foldable / Portable Patient Lift

- For users requiring Maximum Assistance
- 400 lb. (182 kg) safe working load
- Folds flat simply and quickly
- **Locking rear casters**, foot-activated lock/unlock.
- **Floor extraction:** Can be used to lift patients from floor.
- **Quick-release battery** can be charged via charging cradle (optional), or via plug-in charger (included) for charging in place.
- **Overload protection** with visual/audible alert, halts lifting when capacity is exceeded.
- **Manual emergency lowering** for added safety.
- **Simple, fast assembly** using tools provided.



### SPECIFICATIONS

Weight Capacity	400 lb / 182 kg
Overall Base Length	41.0"
External Base Width (Closed)	27.6"
Internal Leg Width (Open)	32.0"
Upright Mast Height	48.3"
Base Height	4.7"
Minimum Lifting Height with Spreader Bar	18.5"
Maximum Lifting Height with Spreader Bar	67.2"
Net Weight	90 lb / 41 kg
Shipping Weight	104 lb / 47 kg

Item #	Description
F400	Foldable and Portable Patient Lift
Warranty	3 years (frame and spreader bar); 2 years (electronics); 1 year (battery)

# Span F350CT

## Car Transfer Patient Lift

- For users requiring Maximum Assistance
- 350 lb. (160 kg) safe working load
- Folds flat for easy transport and storage
- **Locking rear casters**, foot-activated lock/unlock.
- **Floor extraction:** Can be used to lift patients from floor.
- **Quick-release battery** can be charged via charging cradle (optional), or via plug-in charger (included) for charging in place.
- **Overload protection** with visual/audible alert, halts lifting when capacity is exceeded.
- **Manual emergency lowering** for added safety.
- **Simple, fast assembly** using tools provided.



### SPECIFICATIONS

Weight Capacity	350 lb / 160 kg
External Base Width (Closed)	27.6"
Internal Leg Width (Open)	32"
Mast Height	47.3"
Overall Base Length	41"
Minimum Lifting Height	22.2"
Maximum Lifting Height	64.7"
Net Weight	95 lb / 43 kg
Shipping Weight	108 lb / 49 kg

Item #	Description
F350CT	Car Transfer Patient Lift
WWarranty	3 years (frame and spreader bar); 2 years (electronics); 1 year (battery)

