

SPAN For users Requiring SSISTANCE

BEDS Ś SEATING , Mh FALL PROTECTIO

SURFACES

- High Capacity: 500 lb (228 kg) safe working load.
- Care Flexibility: Can be used to lift patient from floor; fits beneath most low beds.
- Safe and Reliable: Redundant controls; overload protection; emergency lowering; 3-year frame warranty.
- Easy operation: Simple charging; audible/visual indicators; locking rear casters.





F500P Full Body Patient Lift

Designed to be the workhorse of your Safe Patient Handling and Mobility program. Delivers a high-value combination of patient safety, caregiver protection, ease of operation and low cost of ownership.





6-point spreader bar maximizes sling options

| SPECIFICATIONS | |
|-------------------------------|-----------------|
| Weight Capacity | 500 lb / 228 kg |
| Base Legs Length | 44.9" |
| Base Closed External Width | 24.9" |
| Base Open External Width | 42.9" |
| Mast Height | 48.3" |
| Base Height | 4.5" |
| Min. lifting height | 25.8" |
| Max. lifting height | 69.1" |
| Net Weight | 99 lb / 45 kg |
| Shipping Weight | 108 lb / 49 kg |



LCD battery level display; audible low-battery indicator



2-button pendant with LED indicator



Underbed clearance, allows use with most low beds

- Easy base width adjustment with locked position; accommodates chairs and commodes.
- **Redundant operation** via control unit, allowing use if pendant is lost or damaged.
- Floor extraction: Can be used to lift patients from floor.
- **Overload protection** with visual/audible alert, halts lifting when capacity is exceeded.
- Manual emergency lowering for added safety.
- Quick-release battery can be charged via charging cradle (optional), or via plug-in charger (included) for charging in place.
- Backlit LCD display during operation for ease of viewing.
- Simple, fast assembly using tools provided.
- Certified to: ISO 10535; IEC/EN/ANSI/AAMI/CAN/CSA-C22.2 60601-1, 60101-1-2, and 60601-1-11; SGS Mark for Safety.

| ltem # | Description |
|-----------|---|
| F500P | Full Body Patient Lift |
| Warranty: | 3 years (frame and spreader bar); 2 years (electronics); 1 year (battery) |

Part of Span's line of Safe Patient Handling and Mobility solutions, including:

F500S Sit-to-Stand Lift

F450T Standing Transfer Aid



Portable Lithium-ion







Overhead Lifts

better science. better care. www.spanamerica.com

800-888-6752 • PO Box 5231 Greenville, SC 29606 PL F500P | CO # 2083 | RV1 • © 2019 Span-America Medical Systems, Inc.

