# PressureGuard® Bariatric APM



Segmented top reduces sharp pressure drop-offs, while air system below handles the "heavy lifting," alternating weight-bearing across the surface in user-selectable intervals. "Two up, one down" inflation pattern ensures maximum patient support throughout the cycle.

### **Powerful Digital Control Unit**

- Up to four times the air output of conventional units.
- Eight electronically-controlled comfort settings.
- Programmable cycle times of 20, 25, 30 or 35 minutes.
- Low pressure indicator and audible alarm.



# **Therapy Provided**

- ✓ Powered Pressure Redistribution
- ✓ Alternating Pressure for users up to 750 lbs.

# **QUICK FACTS**

# Indications:

- Aggressive treatment (through Stage 4)
- Prevention of pressure injuries in immobile or high risk patients
- Can't or won't reposition frequently
- Can be positioned off the pressure injury

# **Key Advantages:**

- Simple operation
- Gentle and effective pressure changes
- Customizable firmness and time settings
- PressureGuard design: stability, safety, anti-shearing

ltem #	Description
AP8042-29	Bariatric APM, 80"L x 42"W x 7"H
AP8048-29	Bariatric APM, 80″L x 48″W x 7″H
AP8054-29	Bariatric APM, 80″L x 54″W x 7″H
* AP8442-29	Bariatric APM, 84″L x 42″W x 7″H
AP8448-29	Bariatric APM, 84″L x 48″W x 7″H
* AP8454-29	Bariatric APM, 84″L x 54″W x 7″H
7500	Bariatric APM Control Unit Only

\* Allow up to 10 working days lead time for this size.

Medicare Approved:	E0277 (Group II)
Weight limit:	750 lbs.
Warranty:	18 months (all components)
Mattress Weight:	Approx. 25 lbs.
UL listed:	Yes

#### Additional bariatric options available – see page 18.

All models are certified for conformance to NFPA 101 (Life Safety Code) – ASTM E1590; Cal Tech Bulletin #117 and16 CFR Parts 1632 and 1633 by an independent testing organization.

