


Span Surfaces for Recessed Deck Beds

Specialty models to accommodate the unique challenges of recessed deck beds.




Featuring Span's proprietary Safety Edge™ two-piece bolster system, adapted to overcome the design restrictions inherent in recessed decks.



Inner bolster "cradle" deflects when patient sits on or moves toward edge of bed . . .



Compressible section to accommodate multiple lengths and foot positions



. . . which helps direct patient's center of gravity safely back toward center of bed.

Choose from among these leading Span surfaces


Geo-Mattress® UltraMax™: Premium therapeutic foam
(Item # UMX85V-29
– see page 24)



PressureGuard® Custom Care®: Non-powered reactive pressure redistribution support surface
(Item # NP85V29
– see page 16)



PressureGuard® Span-Care® Convertible: Non-powered reactive pressure redistribution support surface
(Item # SC85V29
– see page 13)



PressureGuard® Custom Care® Convertible: Convertible treatment with add-on alternating pressure and lateral rotation
(Item # CJ85V29;
add-on control unit
item #6500
– see page 14)



PressureGuard® Custom Care® Convertible LAL: Convertible treatment with add-on microclimate management, alternating pressure and lateral rotation
(Item # CL85V29;
add-on control unit
item #8400
– see page 15)



PressureGuard® Easy Air® LR: Microclimate management, alternating pressure and lateral rotation
(Item # L85VLR-29
– see page 5)

