


# Additional Overlay Options

## Category: Delta Pediatric Overlays

### Delta Foam

Item #	Description	Qty/case
DELTA-C11	Overlay 32"x75"x3.5"	1
DELTA-120	Overlay 34"x75"x3.5"	1
<b>Pediatric sizes:</b>		
DELTA-200	Isolette Matt Overlay 13.5"x25"x1.75"	1
DELTA-202	Warmer Overlay 29"x23.75"x2.25"	1
DELTA-204	Crib Overlay 26.5"x53"x2.75"	1

All Delta foam overlays include Soft Skin.



Outstanding specialty surface for weights below 70 lbs. Ideal for pediatric settings.

### QUICK FACTS

#### Delta Pediatric Overlays

##### Primary Use:

- Skin management in infants and children

##### Key Advantages:


- Proven effectiveness
- High conformance to body contours
- Economical
- Protection in home or hospital

## Category: Preventative Overlays

### Pressure Kair™

Channel cut side-to-side, then sliced head to foot to create individual support cells. Divided into distinctive head, torso, and foot zones.

Item #	Description	Qty/case
SP789-000	Pressure Kair 34"x72"x3.5" with Soft Skin®	4
SP889-000	Pressure Kair 34"x80"x3.5" with Soft Skin®	4



### PRT Matt™

Cross-sliced in head, torso, and foot zones. Distinctive lavender color.

Item #	Description	Qty/case
50960-600	PRT-MATT 34"x72"x3" cross cut with Soft Skin® (fitted)	4
50960-604	PRT-MATT 34"x72"x3" cross cut without Soft Skin®	4




### Span+Guard®

Comfortable, economical overlay. 3" peak, 2" base, in multi-zone design.

Item #	Description	Qty/case
SP799-000	Span+Guard 6 zoned pad, 3"P, 2"B with Soft Skin®	4

(Note: Covered under Medicare code E0199)



### QUICK FACTS

#### Preventative Overlays

##### Primary Use:

- Pressure injury prevention in low-mod risk patients


##### Key Advantages:

- Low cost
- Comfortable
- More conformance to body contours than innerspring mattresses

## Category: Convoluted Overlays

### Convoluted Overlays

Item #	Description	Qty/case
SP22R-000	2" Convoluted Pads, Premium Density	12
SP22S-000	2" Convoluted Pads, Standard Density	12
SP32R-000	3" Convoluted Pads, Premium Density	6
SP35S-000	3" Convoluted Pads, Standard Density	6
SP45R-000	4" Convoluted Pads, Premium Density	4
SP45S-000	4" Convoluted Pads, Standard Density	4



Made from high-quality, combustion modified polyurethane foam. Specifications vary to address any level of surface discomfort. All pads are 72"L x 34"W, and are compression rolled for efficient storage.

### QUICK FACTS

#### Convoluted Overlays

##### Primary Use:

- Comfort

##### Caution:

- Convoluted, or "eggcrate", overlays are recognized for comfort, not pressure management. Call your Span-America representative when selecting the pressure management support surface most effective for your patients' needs.

