

Therapy Provided

- ✓ Microclimate Management (LAL)
- ✓ Alternating Pressure
- ✓ Lateral Rotation¹
- ✓ Powered Pressure Redistribution

QUICK FACTS

Indications

(in addition to those for Easy Air, p. 4):

Lateral Rotation for patients who find frequent manual repositioning painful or need more frequent repositioning

Item #	Description
L7535LR-29	Easy Air LR System 75"L x 35"W x 7"H
L8035LR-29	Easy Air LR System 80"L x 35"W x 7"H
L8435LR-29	Easy Air LR System 84"L x 35"W x 7"H
* L8039LR-29	Easy Air LR System 80"L x 39"W x 7"H
* L8439LR-29	Easy Air LR System 84"L x 39"W x 7"H
* 92726	Easy Air LR System 80"L x 42"W x 7"H
* L8053LR-29	Easy Air LR System 80"L x 53"W x 7"H
* L8453LR-29	Easy Air LR System 84"L x 53"W x 7"H
L85VLR-29	Easy Air LR for recessed deck frame (see page 19)
8200	Easy Air LR control unit only

*Allow up to 10 working days lead time for this size.

Medicare Approved:	E0277 (Group II)
Weight limit:	LR models support 500 lbs in FLOAT and ALTERNATING, 350 lbs in ROTATION.
Warranty:	2 years, (all components)
Mattress Weight:	Approx. 22 lbs.
UL listed:	Yes

Therapy Provided

- ✓ Microclimate Management (LAL)
- ✓ Alternating Pressure
- ✓ Powered Pressure Redistribution

QUICK FACTS

Indications:


(in addition to those for Easy Air, p. 4):

Patients weighing 350-1000 pounds with excess skin moisture (do not use for patients weighing less than 350 lbs.).

Item #	Description
L8042XL-29	Easy Air XL system, 80"L x 42"W x 7"H
L8048XL-29	Easy Air XL system, 80"L x 48"W x 7"H
* L8053XL-29	Easy Air XL System, 80"L x 53"W x 7"H
* L8442XL-29	Easy Air XL System, 84"L x 42"W x 7"H
* L8448XL-29	Easy Air XL System, 84"L x 48"W x 7"H
* L8453XL-29	Easy Air XL System, 84"L x 53"W x 7"H
8300	Easy Air XL Bariatric control unit only

* Allow up to 10 working days lead time for this size.

PressureGuard EASY AIR® LR



The all-in-one Easy Air LR model puts caregiver just a button-push away from a powerful, additional treatment modality, "Lateral Rotation" (LR). The unique PressureGuard approach to LR has been documented to help wounds close more than three times faster than on other commonly provided air therapy surfaces¹. Beyond pure pressure management, Lateral Rotation gives caregivers a versatile tool that is indicated for a wide variety of user needs.

¹Anderson, C. et al. Lateral rotation mattresses for wound healing. *Ostomy/Wound Management* April 2004;50(4):50-62.

WARNING: Cancer and Reproductive Harm – www.P65Warnings.ca.gov

PressureGuard EASY AIR® XL



All models are certified for conformance to NFPA 101 (Life Safety Code) – ASTM E1590; Cal Tech Bulletin #117 and 16 CFR Parts 1632 and 1633 by an independent testing organization.

Additional bariatric options available – see page 20.



SURFACES