



better science. better care.

PressureGuard[®] Air Therapy Surfaces Matrix

Therapy Modalities	Non-powered Air Therapy				Alternating/Lateral Rotation				Microclimate Management ["Low Air Loss" (LAL)]				Bariatric		
	Custom Care [®]	Span-Care [®] Convertible	Custom Care [®] Convertible	Custom Care [®] Convertible LAL	APM	APM ² [™]	Safety Supreme	Span-Care [®] Convertible w/ control unit	Custom Care [®] Convertible w/ control unit	Protocol [®]	Easy Air [®] Easy Air [®] LR	Custom Care [®] Convertible LAL w/ control unit	Bariatric CFT [®]	Bariatric APM	Easy Air [®] XL
Non-powered pressure redistribution	✓	✓	✓	✓									✓		
Convertible to powered via add-on control unit		✓	✓												
Alternating Pressure					✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Lateral Rotation						✓	✓	✓	✓	✓	✓	✓			
Microclimate Management ("LAL")									✓	✓	✓	✓			✓
Powered Flotation					✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Versions to accommodate recessed deck frames	✓	✓	✓	✓											
Alarm: (V)isual, (A)udible					V	V, A	V, A	V, A	V, A	V, A	V	V	V, A	V, A	V
Timed Auto Firm					✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Geo-Matt [®] top surface	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Safety Edge [™] bolster design (RP) = Raised perimeter	✓	✓	✓	✓	✓	✓	✓(RP)	✓	✓	✓	✓	✓	✓	✓	✓
Can maintain air in power outage; (A) = Automatically	✓(A)	✓(A)	✓(A)	✓(A)	✓	✓	✓	✓(A)	✓	✓(A)	✓(A)	✓(A)	✓	✓	✓(A)
Heel Slope [®] , (F) = Foam-only heel section	✓(F)	✓(F)	✓(F)	✓(F)	✓	✓	✓	✓(F)	✓	✓	✓(F)	✓	✓	✓	✓
100% Latex free, 100% DHEP-free	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Non-stretch coated nylon					✓										
LifeSpan [®] bi-directional stretch with Shear Transfer Zone [®] design	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Zippered air flow coverlet with Air Diffusion Matrix [™] layer															
Medicare Approved	E0373	E0373	E0373	E0373	E0277	E0277	E0277	E0277	E0277	E0277	E0277	E0373	E0277	E0277	E0277
Therapeutic Weight Limit (lbs.)	500	500	500	500	350	350	350	400	400	400	500 (350 LR)	750	750	750	1000
Warranty: Mattress, non-pro-rated	5 yrs.	5 yrs.	5 yrs.	5 yrs.	18 mos.	18 mos.	18 mos.	5 yrs.	5 yrs.	18 mos.	2 yrs.	5 yrs.	18 mos.	18 mos.	2 yrs.
Warranty: Control unit, non-pro-rated					18 mos.	18 mos.	18 mos.	2 yrs.	2 yrs.	18 mos.	2 yrs.	2 yrs.	18 mos.	18 mos.	2 yrs.
Meets Federal Fire Safety Regulations	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓