

FGA-700 Dimensions

Name	FGA-700	FGA-700 Low Base
Part Number	280400	280402
Description	FGA-700 Floor Lift	FGA-700 Low Base Floor Lift
Unit	1 piece	1 piece
Lifting Height	25.2-79.4in (640-1940mm)	23.2-74.5in (590-1890mm)
Floor to top of legs	6.1in (155mm)	4.1in (105mm)
Floor to bottom of legs	3.3in (85mm)	1.4in (35mm)
Overall length	59in (1500mm)	
Distance between legs (inside)	25.2-41.9in (640-1065mm)	25.2-41.2in (640-1065mm)
Distance between legs (outside)	30.1-46.8in (765-1190)	30.1-46.1in (765-1170mm)
Min-Max Reach	17.8in-18.9in (453-480mm)	
Max. reach (distance between center of carry bar to front of mast)	27.2in (691mm)	
Turning radius	62.8in (1595mm)	63.1in (1605m)
Dual-wheel castors	Front 4.0in (101mm) Back 5.0in (127mm)	Front: 3.0in (76mm) Back: 4.0in (101mm)
Weight of lift	161lb (73kg)	141lb (64kg)
Safe working load	700lb (318kg)	