

Stand Assist Standard

M921519

- Fully padded with a boned back for added support
- Padded base and under arm rolls for extra comfort
- Hook & loop fastening with clunk-click buckles for patient safety
- Ideal for maintaining mobility and for building strength
- Weight bearing ability required
- Wash Temp 95°C/203°F
- Tumble dry low heat (max 50°C)



Product Codes & Specifications:



100% Polyester		A	B	C
S	M921519-S	87 - 107cm	82cm	22cm
M	M921519-M	97 - 117cm	89cm	22cm
L	M921519-L	112 - 132cm	100cm	23cm
XL	M921519-XL	125 - 145cm	109cm	26cm

