

## Stand Assist - Clips

## M921526

- Fully padded with a boned back for added support
- Padded base and under arm rolls for extra comfort
- Hook & loop fastening with clunk-click buckles for patient safety
- Ideal for maintaining mobility and for building strength
- Weight bearing ability required
- Wash Temp 95°C/203°F
- Tumble dry low heat (max 50°C)

### Product Codes & Specifications:

100% Polyester		A	B	C
<b>S</b>	<b>M921526-S</b>	87 - 107cm	82cm	22cm
<b>M</b>	<b>M921526-M</b>	97 - 117cm	89cm	22cm
<b>L</b>	<b>M921526-L</b>	112 - 132cm	100cm	23cm
<b>XL</b>	<b>M921526-XL</b>	125 - 145cm	109cm	26cm

