⁹ High Easy Sling

WITH LOOPS


- This full body sling transfers in a slightly reclined position and is extremely comfortable and supportive, featuring additional padding on the leg sections.
- The sling is simple and easy to apply and is a go-to sling for many institutions due to its ease of application and minimum options. Includes attached hip tapes as standard for optimum client comfort.
- This sling is suitable for use with the leading standard spreader bar type lifts (loop fixing). Please see the compatibility statement at the back of this guide.
- Available in the following material types: Superfine PLUS[™] (100% polyester spacer fabric; superior comfort over standard sling fabrics; easy application, suitable for bathing.) Mesh (100% polyester open weave mesh; high density, durable, used for bathing as quick drying.)
- * Also available in Poly (M921514-XS; M921514-S; M921514-M; M921514-L; M921514-XL and M921514-XXL).


| // | 7 | | riouuci coucs | |
|---------|------------------------------------|--|--------------------|-------------|
| H | A Shoulder circumference | B Top of head to base of spine | Superfine PLUS™ | Mesh |
| XSmall | 30–34in (76–86cm) | 23.5–26in (60–66cm) | M921521-XS | M921513-XS |
| Small | 34–38in (86–96cm) | 26–29in (66–74cm) | M921521-S | M921513-S |
| Medium | 38–42.5in (96–108cm) | 29–31.5in (74–80cm) | M921521-M | M921513-M |
| Large | 42.5–51.25in (108–130cm) | 31.5–34.5in (80–88cm) | M921521-L | M921513-L |
| XLarge | 51.25–56.75in (130–144cm) | 34.5–39.5in (88–100cm) | M921521-XL | M921513-XL |
| XXLarge | 56.75–63in (144–160cm) | 38–41in (96–104cm) | M921521-XXL | M921513-XXL |



POLY

SUPERFINE PLUS™





Guidelines

• SWL 272kg/600lb

Product Codes

- Wash 95°C/203°F
- Low Tumble Dry