




# In-Situ Deluxe Leg Sling

WITH LOOPS

- The In-Situ Deluxe Leg sling provides full body support on transfer and offers good hip, thigh and femur support. This sling includes integral head support with removable support bones and hip tapes for greater user comfort.
- This divided leg sling is easy to apply, with low-friction fabric to leg areas. The sling includes looped hip tapes offering optimum pick-up options and enhanced user comfort by eliminating pressure from the hip area during lifting.
- In-Situ design may allow sling to be left in place for longer periods with full written risk assessment.
- The sling is perfect for standard transfers but also can be used in bathing transfers as the fabric is a quick drying spacer mesh. Some commode access available, with tailoring to the shoulders, head and commode areas giving a superior fit for even complex body shapes. This sling is suitable for asymmetrical body shapes as the fabric helps to mold to individual body shapes.
- **Superfine PLUS™** (100% polyester spacer fabric; superior comfort over standard sling fabrics; easy application, suitable for bathing.)



**A**




Shoulder circumference

**B**

Top of head to base of spine

**Product Codes**

	<b>A</b>	<b>B</b>	<b>Product Codes</b>
XSmall	30–34in (76–86cm)	23.5–26in (60–66cm)	M921511-XS
Small	34–38in (86–96cm)	26–29in (66–74cm)	M921511-S
Medium	38–42.5in (96–108cm)	29–31.5in (74–80cm)	M921511-M
Large	42.5–51in (108–130cm)	31.5–34.5in (80–88cm)	M921511-L
XLarge	51–56.5in (130–144cm)	34.5–39.5in (88–100cm)	M921511-XL
XXLarge	56.5–63in (144–160cm)	38–41in (96–104cm)	M921511-XXL



## Guidelines

- SWL 272kg / 600lb
- Wash 95°C / 203°F
- Low Tumble Dry