

Universal Sling

M921501

- Full body supportive sling
- Divided padded leg sections
- Ideal for general transfers
- Small aperture for hygiene purposes
- Colour coded tapes
- Wash Temp 95°C/203°F
- Tumble dry low heat (max 50°C)



Product Codes & Specifications:



	100% Polyester Mesh	A	B	C	D
XXXS					
XXS					
XS	M921501-XS	76 - 86cm	44 - 50cm	60cm	54cm
S	M921501-S	86 - 96cm	50 - 57cm	70cm	64cm
M	M921501-M	96 - 108cm	57 - 64cm	90cm	75cm
L	M921501-L	108 - 130cm	64 - 69cm	100cm	81cm
XL	M921501-XL	130 - 144cm	69 - 74cm	110cm	84cm
XXL	M921501-XXL	144 - 160cm	69 - 74cm	130cm	95cm

