



better science. better care.

Seating Matrix Innovative Solutions

	Treatment w/Positioning		Positioning	Skin Protection				General Use					
	Isch-Dish®	Equalizer®	Isch-Dish® Thin	EZ Dish®	Geo-Matt® Contour	Geo-Matt® PRT	Geo-Matt®	Geo-Matt® Wedge	Gel-T®	Soft Pedal	Cross Cut	Short Wave®	Geo-Wave®
General Use					✓			✓	✓	✓	✓	✓	✓
Level of Care	✓	✓		✓	✓	✓							
Positioning Features	✓	✓	✓	✓	✓			✓					
Treatment Surface	✓	✓											
Geo-Matt®					✓	✓	✓	✓				✓	✓
Passive Microclimate Mgmt.					✓	✓	✓	✓				✓	✓
Minimizes Shearing					✓	✓	✓	✓				✓	✓
Gel Components								✓	✓				
Anti-Slip Bottom				✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Machine Washable	✓												
Two-way Stretch Fabric Top	✓	✓		✓									
Wipe Clean				✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Medicare Code *Bariatric Code	E2607	E2607	E2603	E2603	E2603	E2603 E2604*	E2603	E2601	E2601 E2602*	E2601	E2601	E2601 E2611*	
Therapeutic Weight Limit (lbs.)	250	250	250	250	250	250	250	250	250	250	250	300	300
Bariatric Weight Limit (lbs.)		500				500			500				
Warranty	2 yrs.	2 yrs.	2 yrs.	18 mos.	18 mos.	18 mos.	18 mos.	18 mos.	18 mos.	18 mos.	1 yr.	18 mos.	18 mos.
Seat and Back Cushion												✓	✓
Height of Cushion (Inches)	3.5"	3.5"	2.25"	3"	3"	3" standard 4" bariatric	4"	4" front 2" back	2.5"	1.9"	3"	3"	3"
Viscoelastic Top	✓			✓						✓			