


## Bariatric Deluxe, Short Tucks

**M921517**

- Short tucks to back section to allow more fullness to main body of sling
- Elasticated sides
- Full supportive sling
- Manufactured in a woven polyester
- Full leg paddles for femur support with easy glide fabric to both sides for ease of application
- Additional inner leg loops to help facilitate optimum pick up positions
- Access for hygiene purposes
- Wash Temp 95°C/203°F
- Tumble dry low heat (max 50°C)



### Product Codes & Specifications:



100% Heavy Duty Polyester

		A	B	C	D
1	M921517-1	80 - 110cm	up to 75cm	102cm	75cm
2	M921517-2	120 - 140cm	up to 75cm	155cm	75cm

