



Healthcare Facility Design + Project Management Consultancy for Architects and Planners

Your Safe Patient Handling Design and Planning Partner for Hospitals, and Long-Term Care Homes

Handicare has partnered and consulted with architects and builders in the planning and design of hospitals facilities and long-term care homes. Our teams work closely with you to ensure your projects incorporate proper design for all your patient handling needs including ceiling lift and hygiene bathtub installation. It's critical for care teams to have safe patient handling equipment that can be accessed without interference and obstructions, Handicare guides you every step of the way so our products seamlessly fit into spaces for daily care.

We specialize in designing patient lift support systems using the latest Building Information Modelling Services (BIM). Our full-time design and project management teams work directly with our installation technicians to ensure accuracy and quality in the field.

A dedicated Project Manager will be assigned to you to work on all phases of your project, from first assessment, delivery, installation and final acceptance. We have extensive experience working within all health care environments including Long Term Care, Operating and Emergency Rooms, Intensive Care Units, Medical Imaging, Morgues, Therapy, General Medicine and Pool Areas.

We also offer Annual Preventative and Corrective Maintenance on our products through the Handicare Plus program.



Our BIM Services provides:



PE Structural Review

Licensed structural engineer will review attachment and load details and prepare stamped specifications and calculations.



Seismic Design Services

Seismic bracing and vertical reinforcement design and calculations in applicable geographic areas.



MEP Conflict Avoidance

Our BIM teams attend all applicable coordination meetings with MEP trade representation to identify and eliminate clashes.

Custom Solutions for Your Projects

Each new project is unique and presents its own set of challenges. Involving Handicare in the design phase allows you to avoid costly mistakes and meet all your critical timelines. Our ceiling lifts, rail systems, and hygiene bathtubs are speciality pieces of equipment whose functions and operations need to be fully understood before designing room layouts. We take the time to educate you on how our products work, the required space needed, and the best way to implement them into any care scenario, providing you and your clients the best results possible.



Safe Patient Handling Equipment Supplied by Handicare

Handicare Safe Patient Handling Products feature a refined design aesthetic that can complement any place of care they're in. Our products are designed to fit into any room, big or small, and make effective use of all the space that is available to ensure user and patient comfort.

- Ceiling Lifts
- Ceiling Tracks
- Central Bathing Tubs
- Showering
- Floor Lifts
- Beds and Mattresses
- Slings
- Manual Transfer Aids



Change Order Protection

Peace-of-mind knowing that the structure we build has been coordinated to be 100% free of clashes and conflicts.



High Detail Shop Drawings

Install Technicians are prepared, knowing that our shop drawings are accurate down to the fraction of an inch.



As-Built Accuracy

Facility is provided not only the ceiling track locations, but the various structural elements above the ceiling. Great resource for preventative maintenance services.



Proprietary Product Models

Our ceiling tracks and lifts have unique designs and aesthetics. Our custom BIM families accurately represent the appearance and dimensions of our products making a great impression during virtual walkthroughs.

Frequently Asked Questions

When is the right time to include Handicare in the design and planning of our project?

The short answer is right from the start. Safe Patient Handling (SPH) Programs need to be thought of early so they can be planned for your new project. Including our team from the start allows you to optimize your SPH designed spaces with functionality in mind.

What makes up your product portfolio?

Handicare is a full supplier of Safe Patient Handling products that encompass the complete continuum of care. We offer ceiling lifts, track systems, slings, bathtubs, floor lifts, and manual transfer aids. If your project calls for a stairlift or elevator addition, we can assist through our dealer network and parent corporation.


What is your experience working with architects and planners?

We have worked with architects and planners across North America including general contractors and builders. Our teams have also consulted with hospital and nursing home executives, safe patient handling experts, clinical teams, residents, and relatives on the implementation of new and continued projects. Ask your account manager to walk you through some of the unique safe patient handling projects we completed with architects and planners across North America.

HANDICARE

81 Romina Drive, Concord, ON, L4K 4Z9
Toll Free: 877.304.5438 | Fax: 416.260.0256
customerservice.canada@handicare.com
www.handicare.ca

10888 Metro Court, St. Louis, MO 63043
Toll Free: 866.891.6502 | Fax: 314.692.7858
customerservice.northamerica@handicare.com
www.handicareusa.com



Who makes up your project management design team?

The Project Management Team consists of project managers, engineers, construction technicians, and technical drafts staff, all of whom are experts in the implementation of safe patient handling program installations.

Can you provide CAD files for your products?

Yes, we want to make it easy for you to design your care spaces. CAD files are available to use upon request.

Ensure you receive the support you need for your next Safe Patient Handling project. Its never too early to reach out and get started with our Healthcare Facility Design + Project Management Consultancy.

To learn more, call 1-888-637-8155 or visit www.handicare.ca / handicareusa.com