

In-Situ High Easy Clips

M921530

- Designed to provide full body & head support on transfer
- Adjustable integral head support with removable stiffeners
- Easy to apply with low-friction fabric to leg areas for ease of application
- Soft binding as standard
- Universal clip fixing
- Once the patient has been transferred, this sling can be left in-situ, subject to an appropriate risk assessment
- Wash Temp 95°C/203°F
- Tumble dry low heat (max 50°C)



Product Codes & Specifications:



Silva-Superfine PLUS®		A	B	C	D
		66 - 76cm	56 - 62cm	66cm	64cm
XS	M921530-XS	76 - 86cm	60 - 66cm	68cm	69cm
S	M921530-S	86 - 96cm	66 - 74cm	76cm	77cm
M	M921530-M	96 - 108cm	74 - 80cm	86cm	84cm
L	M921530-L	108 - 130cm	80 - 88cm	98cm	93cm
XL	M921530-XL	130 - 144cm	88 - 100cm	102cm	104cm
XXL	M921530-XXL	144 - 160cm	96 - 104cm	120cm	108cm

