

Universal High Back

M921531

- Divided leg sling offering full body & head support
- Full length removable stiffeners
- Manufactured in structured polyester
- Option of a detachable elasticated chest support
- Padding to the leg sections providing comfort
- Wash Temp 85°C/185°F
- Tumble dry low heat (max 50°C)



Product Codes & Specifications:



100% Polyester

		A	B	C	D
XS	M921531-XS	76 - 86cm	60 - 66cm	73cm	81cm
S	M921531-S	86 - 96cm	66 - 74cm	82cm	87cm
M	M921531-M	96 - 108cm	74 - 80cm	92cm	92cm
L	M921531-L	108 - 130cm	80 - 88cm	102cm	95cm
XL	M921531-XL	130 - 144cm	88 - 100cm	117cm	101cm
XXL	M921531-XXL	144 - 160cm	88 - 100cm	130cm	101cm

