

## Universal High Back Mesh

## M921533

- Divided leg sling offering full body & head support
- Full length removable stiffeners
- Manufactured in structured polyester Mesh
- Option of a detachable elasticated chest support
- Padding to the leg sections providing comfort
- Wash Temp 85°C/185°F
- Tumble dry low heat (max 50°C)



### Product Codes & Specifications:



100% Polyester

		A	B	C	D
<b>XS</b>	<b>M921533-XS</b>	76 - 86cm	60 - 66cm	73cm	81cm
<b>S</b>	<b>M921533-S</b>	86 - 96cm	66 - 74cm	82cm	87cm
<b>M</b>	<b>M921533-M</b>	96 - 108cm	74 - 80cm	92cm	92cm
<b>L</b>	<b>M921533-L</b>	108 - 130cm	80 - 88cm	102cm	95cm
<b>XL</b>	<b>M921533-XL</b>	130 - 144cm	88 - 100cm	117cm	101cm
<b>XXL</b>	<b>M921533-XXL</b>	144 - 160cm	88 - 100cm	130cm	101cm

