

## Manual Birthing Sling

- The Manual Birthing Sling is an emergency evacuation sling specifically designed for the birthing pool
- A soft mesh fabric with large holes allows for quick water drainage
- Colour coded grab handles allow for ease of identification and positioning
- This is a manual handling aid and the correct manual handling practices should be followed
- Hoistable version available
- Wash Temp 85°C/185°F
- Tumble dry low heat (max 50°C)

## BIRTH1



### Product Codes & Specifications:



**NB:** 220kg SWL is for the product.  
Strict Manual Handling Regulations must be applied when in use.

	100% Polyester Mesh	A	B	C	D
ONE SIZE	BIRTH1	96 - 104cm	144 - 160cm	110cm	173cm

