


## Leg Lifting Support Sling

- Manufactured in quilted structured fabric
- Offers additional support for patient hygiene & wound treatment
- Helps to reduce patient/carer handling
- Can be used during a transfer or to provide support during a procedure
- Wash Temp 85°C/185°F
- Tumble dry low heat (max 50°C)

## B-LEG



### Product Codes & Specifications:

			
			
	<b>100% Heavy Duty Polyester</b>	<b>A</b>	<b>B</b>
<b>One Size</b>	<b>B-LEG</b>	112cm	21cm

